

Lambda Library User Manual



Clontech

United States/Canada
800.662.2566

Asia Pacific
+1.650.919.7300

Europe
+33.(0)1.3904.6880

Japan
+81.(0)77543.6116

Clontech Laboratories, Inc.
A Takara Bio Company
1290 Terra Bella Ave.
Mountain View, CA 94043
Technical Support (US)
E-mail: tech@clontech.com
www.clontech.com

This protocol has been optimized for use with libraries constructed in the following vectors:

- λ gt10
- λ gt11
- other lambda-derived vectors
- EMBL3
- EMBL3 SP6/T7

Cat. No. 634144
PT1010-1 (PR92374)
Published 18 February 1999

Table of Contents

I. List of Components	3
II. Additional Materials Required	4
III. Library Information	10
A. Library construction	10
B. Library storage	12
C. Titer	12
D. Number of independent clones	12
E. Recombination frequency	12
IV. Bacterial Host Strain Information	13
V. Library Protocols	14
A. Bacterial culture plating	15
B. Library plating/titering	15
C. Probe design	17
D. Library screening using oligonucleotide or DNA probes	18
E. Library screening using an antibody probe	21
F. Isolation of DNA from lambda lysates	26
G. Insert excision	31
VI. Related Products	32
VII. References	33
Appendix A: Host Strain Information	34
Appendix B: Vector Maps	36
Appendix C: Sequence Information	38
1. λ gt10: Nucleic acid sequence around EcoR I cloning site	38
2. λ gt11: Nucleic acid sequence around Sac I–Kpn I region	38
3. λ gt11: Amino acid sequence around Sac I–Kpn I region	42
Appendix D: Recombination Frequency Determination	46
Appendix E: RPM Calculation	47

I. List of Components

Store reagents as described below.

Components included with all λ libraries

- **0.2 ml Liquid cDNA library lysate** (in 1X lambda dilution buffer and 7% DMSO)

- Notes:**
- For long-term storage (6 months or longer), store at -70°C . For short-term storage (less than 6 months), store at 4°C . For additional library storage information, please refer to Section IV.B.
 - All custom libraries are supplied without 7% DMSO. For long-term storage, store phage lysate at -70°C after adding DMSO to 7% or glycerol to 50%; for short-term storage, store at 4°C .

OR

- **0.3 ml Liquid genomic library lysate** (in 1X lambda dilution buffer and 7% DMSO)
- **0.5 ml Bacterial host strain** (in 25% glycerol)

- Notes:**
- All phage libraries are supplied with at least one bacterial host strain. Please refer to Section V for a list of the supplied host strains. For information about bacterial genotypes and references, see Appendix A.
 - Store bacterial cultures (in 25% glycerol) at -70°C . Bacterial cultures stored in this manner should remain viable for more than one year.

II. Additional Materials Required

- **Ampicillin** (50 µg/ml)
Store at -20°C.
- **β-Mercaptoethanol** (14 M)
Store at room temperature (22°C).
- **Carrier DNA**, sheared (10 µg/ml) (λ DNA or calf thymus DNA)
Use of carrier DNA is optional; however, it will increase the transformation efficiency by 5- to 10-fold. Store at -20°C.
- **Cesium chloride** (CsCl)
Store at room temperature (22°C).
- **Chloroform:isoamyl alcohol** (24:1)
Store at room temperature (22°C).
- Biotinylated secondary antibody & avidin-horseradish peroxidase complex or secondary antibody conjugated with horseradish peroxidase & peroxidase substrate solution
- **50X Denhardt's solution**

5.0 g	Ficoll
5.0 g	Polyvinylpyrrolidone
5.0 g	BSA, Pentax Fraction V

Add H₂O to a final volume of 500 ml. Store at -20°C.
- **100% Dimethylsulfoxide** (DMSO)
Store at room temperature (22°C).
- **DNA denaturing solution**

1.5 M	NaCl
0.5 N	NaOH

Store at room temperature (22°C).
- **DNase I**
Store at -20°C.
- **EDTA** (0.5 M)
Store at room temperature (22°C).
- **70% Ethanol**
Store at 4°C.
- **95% Ethanol**
Store at 4°C.
- **Ethidium bromide** (EtBr)
10 mg/ml in H₂O. Store at room temperature (22°C).
- **Fetal calf serum**
Store at 4°C.
- **100% Glycerol**
Store at room temperature (22°C).

II. Additional Materials Required *continued*

- **IPTG** (10 mM & 100 mM)
Isopropyl β -D-thiogalactopyranoside. Filter-sterilize. Store at 4°C.
- **10X Lambda dilution buffer**
58.3 g 1.0 M NaCl
24.65g 0.1 M $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$
350.0 ml 1.0 M Tris-HCl (pH 7.5) (final concentration of 0.35 M)
Add H_2O to a final volume of 1 L. Autoclave and store at 4°C.
- **1X Lambda dilution buffer**
100 ml 10X Lambda dilution buffer
5 ml 2% Gelatin
Add H_2O to a final volume of 1 L. Autoclave and store at 4°C.
The 0.01% gelatin in the 1X lambda dilution buffer stabilizes the library titer for long-term storage. Gelatin is optional when diluting the phage for immediate titering.
- **LB agar plates**
Prepare LB broth as described below. Add agar (15 g/L), and then autoclave. Pour plates, and store them at 4°C.
- **LB agar + 10 mM MgSO_4 plates**
Prepare LB broth as described below. Add MgSO_4 (2.4 g/L) to a 10-mM final concentration and agar (7.5 g/L). Autoclave. Pour plates, and store at 4°C.
- **LB agar + amp + 10 mM MgSO_4 plates**
Prepare LB agar plates as described above. Add MgSO_4 (2.4 g) to a 10-mM final concentration, and then autoclave. After autoclaving, cool LB broth to 50°C before adding ampicillin (50 $\mu\text{g}/\text{ml}$). Pour plates, and store them at 4°C.
- **LB agarose plates**
Prepare LB broth as described below. Add agarose (15 g/L), and then autoclave. Pour plates, and store them at 4°C.
- **LB agarose + 10 mM MgSO_4 plates**
Prepare LB agarose plates as described above. Add MgSO_4 (2.4 g/L) to a final concentration of 10 mM, and then autoclave. Pour plates, and store them at 4°C.
- **LB broth**
10 g/L Bacto-tryptone
5 g/L Bacto-yeast extract
5 g/L NaCl
Adjust pH to 7.0 with 5 N NaOH. Autoclave, and store broth at room temperature (22°C).
- **LB broth + 10 mM MgSO_4**
Prepare LB broth as described above. Add MgSO_4 to a final concentration of 10 mM, and then autoclave. Store broth at room temperature (22°C).

II. Additional Materials Required *continued*

- **LB broth + 10 mM MgSO₄ + 0.2% maltose**

Prepare LB broth + 10 mM MgSO₄ as described above. After autoclaving, cool LB broth to 50°C before adding maltose to a final concentration of 0.2%. Store at room temperature (22°C).

- **LB Soft top agar + 10 mM MgSO₄**

Prepare LB broth as described above. Add MgSO₄ (2.4 g/L) to a 10-mM final concentration and agar (15 g/L), and then autoclave. Store at 4°C.

- **LB Soft top agarose + 10 mM MgSO₄**

Prepare LB broth as described above. Add MgSO₄ to a 10-mM final concentration and agarose (7.2 g/L), and then autoclave. Store at 4°C.

- **MgSO₄ (1 M and 10 mM)**

Store at room temperature (22°C).

- **NaCl (4 M)**

Store at room temperature (22°C).

- **Neutralizing solution**

1.5 M NaCl

0.5 M Tris-HCl (pH 8.0)

Store at room temperature (22°C).

- **Hybridization filters**

Use positively charge-modified, supported nylon filters such as Schleicher & Schnell's Nytran Plus.

- **Nonhomologous salmon sperm blocking DNA**

Shear salmon sperm DNA by passing it 3 times through an 18-gauge needle. Boil the sheared DNA for 10 min. Quickly chill in ice H₂O. Store at -20°C.

- **NucleoBond®, NucleoSpin®, and NucleoTrap® Purification Products**

Clontech offers a complete line of nucleic acid purification products that are useful for purifying library clones and inserts.

- **20%PEG**

20 g PEG 8000 (Fisher Biotech)

90 ml H₂O

Stir gradually at room temperature to dissolve salt. Adjust volume to 100 ml with H₂O. Store at room temperature (22°C).

- **Peroxidase substrate solution**

2 ml 4-chloro-1-naphthol (3 mg/ml in methanol)

10 ml TBS + 0.01 M imidazole

5 µl 30% hydrogen peroxide

Store at 4°C.

- **Phenol:chloroform (1:1)**

Equilibrate with TE buffer (pH 8.0). Store at 4°C, protected from light.

II. Additional Materials Required *continued*

- **Prehybridization solution 1** (for use with DNA probes)

50%	Formamide
5X	SSPE
5X	Denhardt's solution
0.1%	SDS
100 µg/ml	Denatured nonhomologous salmon sperm blocking DNA

 Store at -20°C .
 - **Prehybridization solution 2** (for use with oligonucleotide probes)

6X	SSPE
5X	Denhardt's solution
0.25 %	SDS
100 µg/ml	Denatured nonhomologous salmon sperm blocking DNA
- Note:** Make the prehybridization solutions on the day they will be used. Warm the solutions to dissolve the SDS, and filter through a 0.45–0.8-µm filter. Store at -20°C . You may also use Clontech's ExpressHyb™ hybridization solution (Cat. Nos. 636831, 636832) in place of either prehybridization solution.
- **Proteinase K**
20 mg/ml in H_2O . Filter-sterilize. Store at -20°C .
 - **Radioactive ink**
Mix a small amount of ^{32}P with waterproof black ink.
 - **RNase A** (DNase-free)
Incubate a 10 mg/ml solution of RNase A in [(10 mM Tris-HCl, 15 mM NaCl), pH 7.5] at 90°C for 15 min, and then cool the solution slowly to room temperature. Store in small aliquots at -20°C .
 - **20% SDS**
Store at room temperature (22°C).
 - **Sodium acetate** (3 M)
Store at room temperature (22°C).
 - **20X SSC**

175.3 g	NaCl (3.0 M)
88.2 g	Sodium citrate• $2\text{H}_2\text{O}$ (0.3 M)

 Adjust to pH 7.0 with 10 N NaOH. Add H_2O to a final volume of 1 L. Store at room temperature (22°C).
 - **20X SSPE**

175.3 g	NaCl (3.0 M)
27.6 g	$\text{NaH}_2\text{PO}_4\cdot\text{H}_2\text{O}$ (0.2 M)
40 ml	0.5 M EDTA (final concentration of 0.02 M)

 Adjust to pH 7.4 with 10 N NaOH. Add H_2O to a final volume of 1 L. Store at room temperature (22°C).

II. Additional Materials Required *continued*

- **TBS buffer**
 50 mM Tris-HCl (pH 7.9)
 150 mM NaCl
 Store at room temperature (22°C).
- **TBST**
 TBS + 0.05% Tween-20
 Store at room temperature (22°C).
- **TE buffer**
 10 mM Tris-HCl (pH 7.5)
 1 mM EDTA
 Store at room temperature (22°C).
- **Tetracycline** (15 µg/ml)
 Store at -20°C.
- **Wash buffer 1**
 2X SSC
 0.5% SDS
 Store at room temperature (22°C).
- **Wash buffer 2**
 1X SSC
 0.1% SDS
 Store at room temperature (22°C).
- **Wash buffer 3**
 0.2X SSC
 Store at room temperature (22°C).
- **Whatman 3M paper**
- **X-Gal** (0.1 M)
 Dissolve in DMF. Store at -20°C.

• **Kodak XAR Film**

Additional materials required, by section:

Section IV.B.: Library Storage

- 100% DMSO
- 100% Glycerol

Section VI.A.: Bacterial Culture Plating

- LB broth + 10 mM MgSO₄ + 0.2% maltose
- LB agar plate (without MgSO₄)
- LB broth + 0.2% maltose
- Antibiotics

Section VI.B.: Library Plating/Titering

- LB soft top agar + 10 mM MgSO₄
- LB agar + 10 mM MgSO₄ plates
- 1X Lambda dilution buffer

II. Additional Materials Required *continued*

Section VI.D.: Library Screening Using Oligonucleotide or DNA Probes

- LB broth + 10 mM MgSO₄ + 0.2% maltose
- 1X Lambda dilution buffer
- LB soft top agarose + 10 mM MgSO₄
- Nitrocellulose filters (or nylon membranes)
- Radioactive ink
- DNA denaturing solution
- Neutralizing solution
- Whatman 3M paper
- Prehybridization solutions 1 and 2
- 3X SSC
- Wash buffers 1, 2, and 3
- 20% SDS
- 6X SSPE

Section VI.E.: Library Screening Using an Antibody Probe

- LB agar + 10 mM MgSO₄ plates
- LB agar + amp + 10 mM MgSO₄ plates
- LB broth + 10 mM MgSO₄ + 0.2% maltose
- 1X Lambda dilution buffer
- LB soft top agarose + 10 mM MgSO₄
- 10 mM IPTG
- Nitrocellulose filters
- Radioactive ink
- TBST
- TBST + 20% fetal calf serum
- Biotinylated secondary antibody
- Avidin-horseradish peroxidase complex
- Secondary antibody conjugated with horseradish peroxidase
- Peroxidase substrate solution

Section VI.F.: Isolation of DNA from Lambda Lysates

- LB broth + 10 mM MgSO₄ + 0.2% maltose
- 10 mM MgSO₄
- 1X Lambda dilution buffer
- Chloroform:isoamyl alcohol
- LB agarose + 10 mM MgSO₄ plates
- LB soft top agarose + 10 mM MgSO₄
- LB broth + 10 mM MgSO₄
- DNase I (RNase-free)
- RNase A (DNase-free)
- 20% PEG 8000
- 4 M NaCl
- Cesium chloride
- 0.5 M EDTA
- 20% SDS

II. Additional Materials Required *continued*

- Proteinase K
- Phenol:chloroform
- 3 M Sodium acetate
- 95% Ethanol
- 70% Ethanol
- TE buffer

Section VI.G.: Insert Excision

- NucleoBond[®], NucleoSpin[®], and NucleoTrap[®] Nucleic Acid Purification Products

III. Library Information

A. Library Construction

1. cDNA library preparation

All cDNA libraries are prepared according to a modified Gubler & Hoffman (1983) procedure. 5'-STRETCH PLUS cDNA Libraries are subjected to the most rigorous quality control tests of any commercially available library (see Product Analysis Certificate [PAC] for details).

First- and second-strand cDNA synthesis procedures and post-synthesis handling of the cDNA are optimized to produce libraries with larger inserts, and therefore, a greater probability of full-length transcripts.

5'-STRETCH and 5'-STRETCH PLUS cDNA Libraries were formerly prepared from mRNA that had been completely denatured by methylmercuric hydroxide, to release secondary structures. We have optimized first- and second-strand cDNA synthesis conditions in order to avoid this hazardous compound while maintaining larger insert sizes (1). The procedure uses unique oligo(dT)₂₅d(A/C/G) primers for first-strand cDNA synthesis; the 3'-end nucleotide (A/C/G) positions the primer at the junction of the poly-A tail and encoded transcript. Thus, only a relatively small fraction of the poly-A tail must be reverse transcribed to reach the transcript, and more mRNAs are fully reverse transcribed. Also, specially designed EcoR I hemiphosphorylated adaptors, rather than linkers, are used for cloning (see Section IV.A.4). 5'-STRETCH and 5'-STRETCH PLUS cDNA Libraries prepared using the new procedure are indicated by a catalog number in the 5000 range (e.g., HL5000a, ZL5000b).

After second-strand synthesis, cDNAs with molecular weights less than 400 bp are removed; remaining cDNAs are selected and cloned into lambda vectors.

2. Priming methods

cDNA libraries are oligo(dT)-, random-, or oligo(dT)+random-primed, (see PAC). For oligo(dT)+random-primed libraries, two priming reactions

III. Library Information *continued*

are performed for first-strand synthesis; each reaction is primed with either oligo(dT) or random primers. These reactions are pooled before second-strand synthesis, allowing approximately equal representation of oligo(dT)- and random-primed cDNAs in the same library.

Unidirectionally cloned libraries are primed with an oligo(dT) primer containing a vector-appropriate restriction enzyme site. An adaptor or linker containing the other appropriate enzyme site is ligated to the cDNA before cloning. If an adaptor was used, the cDNA is digested with the enzyme corresponding to the site present on the oligo(dT) primer so that the cDNA can be cloned unidirectionally. If a linker was used, the cDNA is digested with both enzymes before unidirectional cloning.

3. **EcoR I linker/adaptor sequences**

We use an EcoR I linker or specially designed adaptor for libraries containing inserts cloned into the EcoR I site of the vector. If a linker was used, the cDNA is subjected to methylation to protect any internal EcoR I sites. The PAC will state whether or not an adaptor was used.

CCGGAATTCCGG (linker sequence;12-mer)

5'-AATTCGCGGCCGCGTTCGAC-3'

3'-GCGCCGGCGCAGCTG-P-5' (adaptor sequence)

There are advantages to cloning cDNA using the specially designed EcoR I adaptor instead of a linker. No methylation or restriction enzyme digestion of the cDNA prior to cloning is required; therefore, internal EcoR I sequences in the cDNA will not be cut. The EcoR I adaptor also contains flanking Not I and Sal I sites so that the inserts can be removed from the phage using alternative restriction enzymes.

4. **Library amplification**

Unless otherwise stated, all libraries are amplified once. Carefully amplified libraries are faithful copies of unamplified ones. All phage libraries are amplified on solid media in order to avoid unequal growth rates of clones that may occur with liquid media.

If the titer of a library drops below a stable level, the library may be reamplified (to raise the titer) due to the rarity of the particular tissue from which the library was prepared. All reamplified libraries are subjected to the same rigorous quality control as are the once-amplified libraries. The PAC indicates if a library has been reamplified.

5. **Carrier**

In the past, part of Clontech's cDNA library preparation protocol used yeast tRNA as a carrier. Some residual yeast genomic sequences may have been cloned with the cDNA when yeast tRNA was used. The PAC accompanying the library indicates if yeast tRNA was used.

6. **Packaging extract**

McR⁻ packaging extracts are used during library construction.

III. Library Information *continued*

B. Library Storage

1. Short-term

Store working aliquots at 4°C; they should be stable for 6 months.

2. Long-term—premade libraries

We recommend removing 2–3 aliquots (50 µl/aliquot) and storing them at –70°C. The titer of libraries stored at –70°C should remain stable for years. Work from a single tube stored at 4°C to avoid freeze/thaw cycles of the stock sample. You need not remove DMSO prior to use.

3. Long-term—custom libraries

We recommend removing 2–3 aliquots (50 µl/aliquot), adding DMSO to each aliquot to a final concentration of 7%, and storing at –70°C. The titer of libraries stored at –70°C should remain stable for years. Use 1 tube at a time (stored at 4°C) to avoid repeated freeze/thaw cycles of the working stock sample. Alternatively, a sterile solution of glycerol can be added to a final concentration of 50% prior to storage at –70°C. DMSO or glycerol are equally good storage media; neither needs to be removed prior to library use.

C. Titer

The lysate titer is the value obtained when the library was constructed, and it is noted on the PAC. Retiter the library to determine the most accurate titer before screening. For additional information, please see Section VI.B.

Note: If the lysate contains at least 10^8 pfu/ml, the library is representative, and you may proceed with screening procedures. If the titer is below 10^8 pfu/ml, contact our Technical Service Department or your local distributor.

At 4°C, the titer may drop several-fold over a period of one year. Stability will vary between libraries, but low-titer or diluted libraries will be less stable over time. For example, under similar storage conditions, the titer of a library diluted to 10^6 pfu/ml will drop more drastically over time than an equivalent sample with a titer of 10^9 pfu/ml. In some cases, the titer can be stabilized with the addition of gelatin to a final concentration of 0.1 mg/ml.

D. Number of Independent Clones

This number is the number of independent recombinant colonies or independent clones that were in the library before amplification, which is performed to stabilize titer. Most libraries have more than one million independent clones and are representative of the complexity of the cDNA population or genomes involved.

E. Recombination Frequency

Recombination frequencies for λ gt11 libraries are determined after library amplification; for libraries in all other λ vectors, these values are determined

III. Library Information *continued*

before amplification. Although the recombination frequency value that we obtained is noted on the PAC, we recommend that you replat the library to determine this value most accurately. See Appendix D for recombination frequency determination for expression libraries—i.e., libraries with selection criteria such as blue/white selection in λ gt11.

IV. Bacterial Host Strain Information

Each library is supplied with at least one bacterial host strain for library plating. Store host strains at -20°C when not in use. For information on bacterial host strain genotypes, please refer to Appendix A.

- Notes:**
- Antibiotics are required for stock maintenance only and should **not** be used with immunoscreening protocols.
 - Use LB agar plates that do not contain MgSO_4 for stock maintenance.

TABLE I: HOST STRAIN / ADDITIONAL INFORMATION

Vector	Bacterial Host	Additional Information
λ gt10	C600 Hfl C600	3, 6 4
λ gt11	Y1090r ⁻	1, 2, 5
EMBL3	K803	
EMBL3 SP6/T7	K802	
EMBL4	K803	

Additional Information:

1. K802 and LE392 may also be used for titering, but not for screening.
2. This vector will form tiny plaques in Y1090r⁻.
3. This strain is recommended for recombinant screening.
4. This strain is recommended for vector propagation. It will also show nonrecombinants in unamplified libraries.
5. Add 50 $\mu\text{g}/\text{ml}$ ampicillin.
6. Add 15 $\mu\text{g}/\text{ml}$ tetracycline.

V. Library Protocols




Figure 1. Overview of the Lambda Library Protocols.

V. Library Protocols *continued*

A. Bacterial Culture Plating

1. Streak 5 μ l of the appropriate E. coli host strain (Table I) from the 25% glycerol stock culture provided onto a MgSO_4 -free LB agar plate containing the appropriate media additives (Table I). Incubate at 37°C overnight. This plate is the master plate and should be stored at 4°C. Fresh streak plates should be prepared at 2-wk intervals.
2. Pick a single, isolated colony, and streak onto another MgSO_4 -free LB agar plate containing the appropriate media additives. Incubate at 37°C overnight. This plate is the primary working plate and should be stored at 4°C.

B. Library Plating/Titering

1. Pick a single, isolated colony from the primary working plate (Step VI.A.2.), and inoculate LB broth + 10 mM MgSO_4 + 0.2% maltose (without antibiotics). Incubate on a shaker (200 rpm) at 37°C overnight until the OD_{600} of the culture reaches 2.0.
2. Preparation of titering dilutions of phage lysate (library)
 - a. Pipet 2 μ l of the library lysate into 1 ml of 1X lambda dilution buffer (Dilution 1 = 1:500).
 - b. Transfer 2 μ l of Dilution 1 into a second tube containing 1 ml of 1X lambda dilution buffer (Dilution 2 = 1:250,000).
 - c. Prepare 4 tubes as described in Table II using the bacterial overnight culture obtained from Step B.1.

TABLE II: LIBRARY PLATING DILUTIONS

Tube	1X Lambda Dilution Buffer	Bacterial Overnight Culture	Phage Dilution 2
1	100 μ l	200 μ l	2 μ l
2	100 μ l	200 μ l	5 μ l
3	100 μ l	200 μ l	10 μ l
4 (Control)	100 μ l	200 μ l	0 μ l

3. Incubate tubes in a 37°C water bath for 15 min.
4. Add 3 ml of melted LB soft top agar + MgSO_4 (45°C) to each of the 4 tubes. Mix well.

V. Library Protocols *continued*




Figure 2. Library screening using a DNA or oligonucleotide probe.

V. Library Protocols *continued*

5. Pour the contents from each tube onto 4 separate 90-mm LB agar + MgSO₄ plates. Swirl the plates quickly after pouring to allow even spreading of the agar.

Note: If the plates were prepared fresh and have been sitting at room temperature, prewarm plates at 37°C for 10–15 min. If plates have been stored at 4°C, prewarm plates at 37°C for 1 hr prior to use.

6. Cool plates at room temperature for 10 min to allow the inoculum to soak into the agar.
7. Incubate plates (inverted position) at 37°C for at least 6–7 hr.
8. Count the plaques to determine titer (pfu/ml). This plate is the secondary working plate and should be stored at 4°C.

$$\text{pfu/ml} = \frac{\text{\# of plaques}}{\mu\text{l used}} \times \text{dilution factor} \times 10^3 \mu\text{l/ml}$$

C. Hybridization Screening—General Considerations

Lambda libraries can be screened by hybridization using either DNA probes or oligonucleotide probes (Section VI.D). In either case, we recommend that you prepare at least two replicate filters of each plate to be screened (Section VI.C). Along with your library screening filters, perform a complete hybridization experiment with a blank filter (or preferably a control filter containing nonrecombinant vector plaques) to determine the conditions under which background noise is lowest. The protocol given below is for screening with a radioactively labeled probe and detecting positive signals by autoradiography.

Further information on hybridization using nucleic acid probes can be found in the laboratory handbooks edited by Sambrook et al. (2001) and Ausubel et al. (1987-1995). Also see Jacobs et al. (1988), Lathe (1985), Suggs et al. (1981), Wallace et al. (1981), and Wood et al. (1985) for further information on using oligonucleotide probes, and Woo (1979) for further information on using DNA probes.

Important notes:

- Use soft top **agarose** (rather than soft top agar) in the following procedures because commercially supplied agar contains components that may interfere with hybridizations.
- Always swirl plates quickly after pouring melted top agarose to promote even spreading of the agarose over the plate.
- Use positively charge-modified, supported nylon filters (Schleicher & Schuell's Nytran Plus works well). Nitrocellulose filters are too brittle for use in hybridization screening.
- The hybridization probe (DNA or oligonucleotide) should contain at least 10⁷ cpm. The specific activity of the probe should be >8 x 10⁸ cpm/μg, and the final concentration 1–2 x 10⁶ cpm/ml of hybridization solution.

V. Library Protocols *continued*

- Commercially available rapid hybridization solutions (e.g., ExpressHyb Hybridization Solution, (Cat. No. 636831) may be substituted for standard prehybridization and hybridization solutions. If you are using a rapid hybridization solution, follow the protocol that is sent with that product.

λ gt11 libraries

Y1090r⁻, the host strain for λ gt11 libraries, contains copies of pMC9, the parental plasmid of pBR322. In amplified λ gt11 libraries, there is a 10^{-4} to 10^{-5} probability that the pMC9 sequence will recombine with the λ gt11 vector. Therefore, when you are screening an amplified λ gt11 library, we recommend that your probe contain no pBR322 or pBR322-like sequences.

Studies have indicated that some λ gt11 inserts may undergo homologous recombination with the vector arms (9). Therefore, false positives may result. We recommend that all positive clones be verified after identification.

D. Library Screening Using Oligonucleotide or DNA Probes

Note: It is essential to use agarose in the following procedures rather than agar as the top layer; agar may peel off the plate during the filter lift. Always swirl plates quickly after pouring to allow for even spreading of the agarose over the plate.

1. Filter replicas

- a. Pick a single, isolated bacterial colony from the primary working plate prepared in Step VI.A.2., and inoculate 15 ml of LB broth + 10 mM MgSO_4 + 0.2% maltose. Incubate with good aeration at 37°C overnight (or until stationary phase is reached).
- b. Based on the titer obtained in Section VI.B., prepare a sample of the phage library diluted in 0.1 ml of sterile 1X lambda dilution buffer so that the dilution will yield a maximum of approximately 10,000 pfu for each 90-mm plate (or approximately 30,000 pfu for each 150-mm plate). Combine the dilution with 200 μ l of bacteria grown in Step a above (or 600 μ l for a 150-mm plate). Incubate at 37°C for 15 min.
- c. Immediately add 3 ml of melted LB soft top agarose + 10 mM MgSO_4 to the cell suspension for each 90-mm plate (or 7 ml for each 150-mm plate). Mix and pour onto an LB agar + MgSO_4 plate prewarmed in a 37°C incubator.

Note: The LB agar plates must be dry, otherwise the plaques may run together. Use 2-day-old plates that have been prewarmed to 37°C for at least 4–5 hr.

- d. Invert the plates and incubate at 37°C until the plaques reach a diameter not exceeding 1.5 mm or are just beginning to make contact with one another (3–8 hr). Avoid confluent lysis.

Note: Plaque densities vary among λ phage and their plaque morphologies. Optimize the number of plaques per plate according to your specific phage.

V. Library Protocols *continued*

- e. Chill the plates at 4°C for at least 1 hr to allow the LB soft top agarose to harden. Plates can now be stored at 4°C overnight.
 - f. Number a nylon membrane with a soft pencil or ball point pen. Using sterile forceps, place the filter onto the LB soft top agarose. Avoid trapping air bubbles. Mark the filter in 3 asymmetric locations by stabbing through the filter and into the agar with an 18-gauge needle attached to a syringe containing radioactive ink.
 - g. After 2 min, carefully peel off the filter. Float the filter on top of DNA denaturing solution, plaque-side up, for 30 sec. Then immerse the filter for 5 min.
 - h. Remove the filter and immerse in neutralizing solution for 5 min.
 - i. Briefly rinse the filter in 2X SSC, and place on a Whatman 3M paper to dry.
 - j. Place a second membrane onto the same plate, and mark it with ink at the same locations. Peel the filter off after 3 min. Denature the DNA; neutralize and rinse as described in Steps g–i. Additional membranes can be placed on the same plate and peeled off after 5 min.
 - k. Use a UV crosslinker to fix the DNA to the filter. Follow the manufacturer's recommended conditions, or simply bake membranes at 80°C for 2 hr.
 - l. Proceed to Step 2 for hybridization using DNA probes or to step 3 for hybridization with oligonucleotide probes.
2. Hybridization using DNA Probes (>200 bp)
- a. Incubate the filters in prehybridization solution I in a heat-sealed plastic bag at 42°C for 4 hr with agitation. Avoid trapping air bubbles in the bag. The filters should come in direct contact with the solution.
Note: Commercially available rapid hybridization solutions (e.g., ExpressHyb™ Hybridization Solution, Cat. No. 636831) may be substituted for standard prehybridization and hybridization solutions.
 - b. Denature the labeled DNA probe by heating at 100°C for 10 min. Chill it quickly on ice, and then add it to the prehybridization solution in the plastic bag. Incubate the membranes in at least 10 ml of hybridization solution (prehybridization solution I containing probe) at 42°C for 16–20 hr (13). If using ExpressHyb™ solution, follow the recommendations provided for library screening.
Note: The probe should contain at least 10⁷ cpm.
 - c. After hybridization, carefully and safely remove the filter from the plastic bag. Wash the filters at room temperature in wash buffer 1 for 15–20 min. Do not allow the filters to dry.

V. Library Protocols *continued*

- d. Wash the filters in wash buffer 2 at 65°C for 1 hr. Repeat this wash once. If the background is still too high, you may wish to increase the washing stringency by incubating the filters for 1 hr in wash buffer 3.

Note: We recommend performing a complete experiment with a blank filter to determine the conditions under which background noise is lowest.

- e. Wrap the filters in plastic wrap, and place them on Whatman 3M paper at room temperature. Pair duplicate filters, and place them onto a piece of filter paper backing cut to the size of the x-ray film to be used for autoradiography.
- f. Mark the filter paper backing at 3 asymmetric locations with radioactive ink. Cover the filters in plastic wrap. Proceed with autoradiography at -70°C for 16–24 hr, using an intensifying screen to enhance the signal. We recommend using Kodak XAR Film.

Note: If the filters are kept slightly moist in the plastic wrap, they can be washed again to reduce background noise.

- g. After developing, align film with the filters to find positive plaques. In cases where a single, isolated positive plaque cannot be picked, remove an agar plug containing several plaques into 1 ml of sterile 1X lambda dilution buffer. Replate to obtain 200–1,000 plaques on a 150-mm plate. Rescreen these plaques. Pick a single, well-isolated plaque for the plate stock. The plate stock is then used for further analysis. For information on the isolation of DNA from lambda lysates, please refer to Section VI.F.

3. Hybridization using oligonucleotide probes (6, 7, 11, 12, 14)

- a. Incubate the filters in prehybridization solution 2 in a heat-sealed plastic bag at 20°C below the oligonucleotide probe's estimated dissociation temperature (T_d) for at least 4–6 hr. Avoid trapping air bubbles in the bag; the filters should come in direct contact with the prehybridization solution. For probes 14–27 nucleotides in length, the approximate T_d is:

$$T_d (^{\circ}\text{C}) = 4(\text{G}+\text{C}) + 2(\text{A}+\text{T})$$

under standard conditions of 1.0 M Na⁺ (6X SSPE or 6X SSC)

For longer probes:

$$T_d = 81.5^{\circ}\text{C} + 16.6 \log_{10} [\text{Na}^+] + 0.41(\% \text{G/C content}) - 500/n$$

n = the number of nucleotides

The T_d of the oligonucleotide can be decreased by approximately 13°C when the Na⁺ concentration of the incubation solution is decreased from 1 M (6X SSPE) to 0.2 M (1X SSPE). 1X SSPE in the prehybridization and hybridization solutions should be used when employing oligonucleotides with a T_d higher than 80°C.

V. Library Protocols *continued*

- b. Denature the oligonucleotide probe by heating at 100°C for 2 min (to release any secondary structure which may have formed). Chill probe quickly on ice. Add the denatured probe to the prehybridization solution in the plastic bag. Incubate the filters in hybridization solution (prehybridization solution containing probe) at the appropriate temperature for 16–20 hr.
- c. After hybridization, carefully discard the solution, and wash the filters under low stringency conditions in a large volume of 2X SSC and 0.05% SDS for 1 hr. The wash buffer should be at room temperature or 30°C below the calculated T_d , whichever is higher. Replace the wash solution with fresh solution at a temperature 20°C below the calculated T_d . Wash the filters under these high stringency conditions for 15 min. If the background is still too high, increase the high stringency wash time to 30 min.

Note: We recommend performing a complete experiment with a blank filter to determine the conditions under which background noise is lowest.

- d. Dry the filters by blotting on Whatman 3M paper at room temperature. Pair duplicate filters, and tape the dried filters onto a piece of filter paper backing, cut to the size of the x-ray film to be used for autoradiography. We recommend using Kodak XAR Film.
- e. Mark the filter paper backing at 3 asymmetric locations with radioactive ink. Cover the filters in plastic wrap. Proceed with autoradiography at –70°C for 16–24 hr, using an intensifying screen to enhance the signal.
- f. After developing, align film with the filters to select a positive plaque. In cases where a single, isolated positive plaque cannot be picked, remove an agar plug containing several plaques into 1 ml of sterile 1X lambda dilution buffer. Replate to obtain 200–1,000 plaques on a 150-mm plate. Rescreen these plaques. Pick a single, well-isolated plaque for the plate stock. The plate stock is then used for further analysis. For information on the isolation of DNA from lambda lysates, please refer to Section VI.F.

E. Library screening using an antibody probe

(Figure 2; Refs 2, 5, 16–18)

Cloning directly into the lacZ gene of λ gt11 results in the expression of foreign DNA as part of a β -galactosidase fusion protein. λ gt11 expresses the α -subunit of β -gal (96 kd). After plating, proteins released by the lysed cells are immobilized onto a nitrocellulose filter overlay. Your primary antibodies (rabbit, goat, or other) are used as probes to detect recombinants containing DNA sequences of interest by probing for specific antigens. Primary antibody binding is visualized by the use of a secondary antibody followed by colorimetric or other detection.

V. Library Protocols *continued*




Figure 3. Library screening using an antibody probe.

V. Library Protocols *continued*

This protocol describes an avidin-horseradish peroxidase (avidin-HRP)-based system for colorimetric detection of the protein products expressed by the cloned cDNA. Antibody dilutions and other reagent must be optimized empirically or according to the manufacturer's recommendations.

1. Pick a single, isolated colony from the primary working plate (Step VI.A.2.), and grow the culture to saturation in 15 ml of LB broth + 10 mM MgSO_4 + 0.2% maltose at 37°C with good aeration.

Note: Maltose addition is required for proper receptivity of cells to phage.

2. Mix 0.2 ml of the Y1090r⁻ culture with 0.1 ml of sterile 1X lambda dilution buffer containing 10^4 pfu of λ gt11 phage. Allow adsorption of phage to the cells at 37°C for 15 min.

Note: Empirically determine the number of λ gt11 clones to plate. 10^4 plaques/90-mm plate is generally acceptable. More than 10^4 plaques may cause a gradual decrease in the intensity of positive signals.

3. Add 3 ml of LB soft top agarose + 10 mM MgSO_4 for each 90-mm plate (or 7 ml for each 150-mm plate). Mix and pour onto an LB agar + 10 mM MgSO_4 plate. Swirl quickly while pouring to allow thorough mixing. Incubate plates at 42°C for 3.5 hr to allow lytic phage growth.

Note: Do not use agar for the top layer of the plate.

4. Using forceps, place a dry nitrocellulose filter previously saturated in 10 mM IPTG on the top agarose. Avoid trapping air bubbles. Mark the filter in 3 asymmetric locations by stabbing through the filter and into the agarose with an 18-gauge needle attached to a syringe containing waterproof ink.

Notes:

- Do not substitute nylon filters for the nitrocellulose filters.
- Filters must be dry so that moisture does not smear the plate. Air-drying the filters at room temperature is best.
- IPTG, a gratuitous inducer, is used to enhance expression of the β -galactosidase fusion gene at 37°C. In the case of the chicken ovalbumin clone (cl), a 3-fold stimulation is detected.

5. Incubate the plates containing filters at 37°C for 3.5 hr.

Notes:

- Do not incubate overnight because increased background may result.
- If a duplicate is required, simply overlay a second filter on top of the agarose after the first filter is removed. Incubate the second filter for an additional 3.5 hr.

6. Remove and rinse the filters in TBST.

Note: If agarose sticks to the filter, cool the plates at 4°C for 15 min before lifting the filters.

V. Library Protocols *continued*

7. Incubate the filters in TBST + 20% fetal calf serum for 30 min. Use 5 ml for each 82-mm filter, or 10 ml for each 132-mm filter.

Note: Fetal calf serum is a blocking agent used to saturate nonspecific protein binding sites on the nitrocellulose filter and reduce background.

8. Remove and rinse the filters in TBST.

NOTE: YOU MAY STOP HERE. STORE COMPONENTS AT 4°C.

9. Primary Antibody binding

- a. Incubate the filters in 5 ml for each 82-mm filter (or 10 ml for each 132-mm filter) of TBST + the appropriate dilution of your primary antibody at room temperature for 1 hr.

- b. Add 2 μ l of your primary antibody/10 ml of solution used above.

Note: The amount of primary antibody may have to be optimized.

10. Wash the filters in 3 changes of TBST; allow 3 min/change.

11. Detection of bound primary antibody with a biotinylated secondary antibody: An alternative detection procedure is described in Step 12. Using a biotinylated antibody will result in stronger signals, albeit with higher background.

- a. Add 20 μ l of secondary antibody/10 ml of TBST. Incubate at room temperature for 30 min with gentle agitation.

- b. Place filters in the TBST containing biotinylated secondary antibody.

- c. Wash the filters in 3 changes of TBST; allow 3 min/change.

- d. Transfer the filters to TBST containing avidin-horseradish peroxidase complex. Incubate at room temperature for 30 min with gentle agitation.

Note: Prepare avidin-horseradish peroxidase (HRP) complex by adding avidin and horseradish peroxidase conjugate to 10 ml of TBST. Consult the manufacturer's recommendations for amounts of these reagents to add. Incubate this mixture at room temperature for 30 min before use.

- e. Wash the filters 3 times in TBS (without Tween-20); allow 3 min/change.

- f. Incubate the nitrocellulose filters in peroxidase substrate solution. Allow the color to develop for approximately 30 min. Use 5 ml of peroxidase substrate solution for each 82-mm filter (or 10 ml for each 132-mm filter). Consult the manufacturer's recommendations to determine the optimal concentration of peroxidase substrate.

Note: Positive plaques appear blue as a result of horseradish peroxidase staining. Once you locate them, remove them from the original plates, and purify recombinants by repeating the same screening procedure until all plaques are positive.

V. Library Protocols *continued*

- g. Wash the filters in three changes of distilled H₂O, and allow the filters to air-dry.

Notes:

- The primary antibody is bound specifically by the biotinylated secondary antibody whose biotin moieties can be easily detected by avidin-conjugated horseradish peroxidase followed by immunostaining.
- The peroxidase stain on the filter is light sensitive. Store the filter in the dark to prevent fading.

12. Detection of bound primary antibody: Use of goat anti-rabbit antibody conjugated with horseradish peroxidase

Note: This detection procedure is alternative to that described in Step 11. You will obtain less background and slightly weaker positive signals when using goat anti-rabbit antibody conjugated with horseradish peroxidase for detection of primary antibody.

- a. Add 8 μ l of affinity-purified horseradish peroxidase-conjugated goat anti-rabbit IgG/10 ml of TBST.
 - b. Transfer the filters to the TBST containing IgG.
 - c. Incubate at room temperature for 1 hr with gentle agitation.
 - d. Wash filters 3 times in TBS (without Tween 20); allow 3 min/change.
 - e. Incubate the nitrocellulose filters in peroxidase substrate solution. Allow color development for approximately 30 min. Use 5 ml of peroxidase substrate solution for each 82-mm filter (or 10 ml for each 132-mm filter).
 - f. Wash with three changes of distilled H₂O, and allow the filters to air-dry.
13. Align the filters with the plates to select a positive plaque. If you cannot pick a single, isolated positive plaque, remove an agar plug containing several plaques, and place it in 1 ml of sterile 1X lambda dilution buffer. Replate to obtain 200–1,000 plaques on a 150-mm plate and rescreen. Pick a single, well-isolated plaque for the plate stock, and use this stock for further analysis. For information on the isolation of DNA from lambda lysates, please see Section VI.F.

V. Library Protocols *continued*

F. Isolation of DNA from Lambda Lysates

1. Preparation of bacterial host cells

- a. Incubate the appropriate bacterial host cells in 15 ml of LB broth (pH 7.5) + 10 mM MgSO_4 + 0.2 % maltose at 37°C overnight.
- b. Pellet the cells, and resuspend the pellet in 7.5 ml of 10 mM MgSO_4 . Store the cells at 4°C for a maximum of 2–3 days following their preparation.

2. Phage stock preparation

Note: The ratio of phage to bacterial host, which determines the ultimate success of any phage DNA preparation, must be determined empirically (9). A yield of 100–200 μg of phage DNA can be expected from 500 ml of phage lysate.

- a. Plate an appropriate number of lytic phage (usually 600–800 plaques per 90-mm plate) from the plate stock from Section VI.D.2.g., D.3.f., or E.13., so that a single plaque can be easily removed. Incubate at 37°C overnight.
- b. With the end of a pasteur pipette, transfer an agar plug containing a single plaque to a microcentrifuge tube containing 200 μl of 1X lambda dilution buffer.
- c. Add a drop of chloroform and vortex briefly.
- d. Let phage elute at 4°C overnight (or at 37°C with shaking [200 rpm] for 4–6 hr).

Note: A low titer may result if phage is allowed to elute at 4°C for longer than 12 hr.

- e. Spin in a microcentrifuge at 10,000 rpm (8,000 x g) for 2 min to remove debris.
- f. Titer the supernatant as described in Section VI.B. Plate 10^5 pfu on a 150-mm LB agarose + 10 mM MgSO_4 plate.
- Note:** Do not use agar for the plate as commercially available agar often contains strong inhibitors of restriction enzymes.
- g. Based on the titer obtained above, plate enough phage so that near confluency is obtained in 5–7 hr. Incubate at 37°C for 5–7 hr.
- h. Add 10 ml of 1X lambda dilution buffer to the plate, and incubate at 4°C overnight.
- i. Add a few drops of chloroform to the plate and swirl briefly. Pour the liquid (containing phage lysate) from the plate into a sterile 50-ml polypropylene tube.
- j. Add 2 ml of chloroform to the plate lysate, and vortex for 2 min.
- k. Centrifuge in a Beckman J2-21 centrifuge at 7,000 rpm (7,200 x g) for 10 min. Collect and save the supernatant.

Note: This is a high-titer stock. Check the titer for future reference as described in Section VI.B. Titer should be approximately 10^{10} pfu/ml.

V. Library Protocols *continued*

3. Lysate preparation

Note: For preparation of large quantities of λ DNA, the liquid lysate method is more convenient than the plate lysate method; however λ DNA generally grows better in plate lysates.

a. Plate lysate

Note: Agarose is recommended for both the top and bottom layers of the plate since commercially available agar often contains strong inhibitors of restriction enzymes.

The number of plates required depends upon intended usage (i.e., preparations of 0.5–1 L of lysate requires 50–100 150-mm plates).

i. Based on the titer of the high-titer stock, prepare lysate dilutions such that approximately 10^5 pfu/ml are obtained on each 150-mm LB agarose + 10 mM MgSO_4 plate.

(a) Pick a single, isolated colony from the primary working plate (Step VI.A.2.), and inoculate LB broth + 10 MgSO_4 + 0.2% maltose. Incubate on a shaker (200 rpm) at 37°C overnight.

(b) Set up the appropriate number of tubes with 600 μl of bacterial culture and diluted lysate.

(c) Incubate in a 37°C water bath for 15 min.

(d) Add 5.5–6.5 ml of melted LB soft top agarose + 10 mM MgSO_4 .

ii. Mix and pour inoculum onto 150-mm plates. Swirl the plates quickly while pouring to allow even spreading of the agarose.

iii. Incubate plates at 37°C for 5–6 hr.

Note: If the transfection is performed in the afternoon, incubate the plates overnight. However, daytime incubation is preferred so that the plates can be harvested at the first sign of confluency.

iv. Add 12 ml of 1X lambda dilution buffer to each plate. Store plates at 4°C overnight.

v. Incubate at room temperature for 1 hr with constant shaking.

vi. Remove and save the 1X lambda dilution buffer solution; rinse the plate surface with an additional 2 ml of 1X lambda dilution buffer.

vii. Pool the lambda dilution buffer solutions. This is the plate lysate.

viii. For lysate processing directions, proceed to Section VI.F.4.

b. Large-scale liquid lysate

Most liquid lysate protocols start by adding a small number of infected cells to a large volume of culture media (14). We have found that better results are achieved with a scaled-up version beginning with a higher quantity of infected cells.

i. Add 1–3 ml of phage stock ($1\text{--}3 \times 10^{10}$ pfu) to 1 L of host cells grown in LB broth to an OD_{600} of 0.6.

Note: Calibrate the spectrophotometer by taring with an LB broth blank.

V. Library Protocols *continued*

- ii. Shake the culture at 37°C in a 4-L flask until lysis is apparent (6–10 hr, depending upon the vector).

Notes:

- When lysis occurs, the culture looks clear with some bacterial debris. If lysis does not occur, the culture appears turbid.

- Foaming during this period is a good indication of phage growth.

- iii. After lysis, add 10 ml of chloroform, and shake for an additional 15 min.

Note: Chloroform will lyse cells which were not lysed during the previous incubation.

- iv. Remove the flask. Centrifuge the lysate at 7,000 rpm (7,200 x g) for 10 min at 4°C in polypropylene bottles.

- v. Combine the supernatants, and check the titer as in Section VI.B. This is the liquid lysate which is now ready for DNA extraction.

Note: If the titer is below 10^8 pfu/ml, repeat Steps i–v. For large quantities of DNA, the titer must be greater than 10^{13} pfu/L.

Standard yield = 440 µg of DNA/L:

$$(10^{13} \text{ pfu/L}) \times (40 \text{ kb/pfu}) \times (660,000 \text{ g/6} \times 10^{23} \text{ kb}) = 440 \text{ µg of DNA/L}$$

If the titer is insufficient, increase the ratio of phage to bacteria.

If lysis appears complete within 2–3 hr, the ratio of phage to bacteria is too high, and the bacterial population has been lysed prematurely. If lysis occurs too quickly, repeat Steps 3.b.i–v.

4. Lysate processing

- a. Centrifuge the bottled, pooled lambda dilution buffer solutions (obtained either from plate or liquid lysates) in a Beckman J2-21 at 8,500 rpm (10,000 x g) for 10 min in order to pellet debris.
- b. To the supernatant, add DNase I to 1 µg/ml and RNase A to 5 µg/ml. Incubate at room temperature for 30 min.
- c. Add 100% chloroform to a final concentration of 5%, and vortex for 30 sec.
- d. Centrifuge at 8,500 rpm (10,000 x g) for 10 min at 4°C to pellet the debris. Transfer the aqueous phase to a new centrifuge tube.
- e. Add an equal volume of 20% PEG/2.0 M NaCl to the aqueous phase. Incubate on ice for at least 1 hr.

Note: The solution may be stored at 4°C overnight.

- f. Centrifuge the precipitated phage at 8,500 rpm (10,000 x g) for 15 min at 4°C. Discard the supernatant. A grayish smear should be evident on the side of each bottle.

Notes:

- Remove as much of the PEG solution as possible as PEG can inhibit restriction enzymes.

- For rapid DNA preps, Steps g–l may be omitted. However, ultracentrifugation is recommended before preparative digestion of cDNA inserts for making probes or subcloning.

V. Library Protocols *continued*

- g. Resuspend the pellets in 32 ml of 1X lambda dilution buffer.
- h. Transfer the phage suspension into two 50-ml polypropylene tubes, and add an equal volume of chloroform. Vortex for 30 sec.
- i. Centrifuge in a Beckman J2-21 centrifuge at 7,000 rpm (7,200 x g) for 10 min. Collect the supernatant, being careful to leave the PEG interface behind.
- j. Add 0.5 g of CsCl per ml of phage suspension.
- k. Pour the phage suspension into a 40-ml UltraClear™ tube, and centrifuge in a Beckman SW 28 rotor at 22,000 rpm (90,000 x g) for 2 hr at 20°C.

Note: The phage will pellet to the bottom of the tube.

- l. Pour off the supernatant, and resuspend the clear, sticky phage pellet in 1 ml of 1X lambda dilution buffer.
- m. Transfer to a 1-ml microcentrifuge tube, and spin down debris at 15,000 rpm for 10 min.

Note: Unless purity is critical (e.g., for making vector arms), the phage is now ready for DNA extraction (proceed to Step 5 below). If purity is essential, proceed to Step n.

5. DNA extraction

Note: You may wish to dialyze rather than precipitate the DNA because large DNAs are often difficult to resuspend once precipitated.

- a. To the phage DNA, add EDTA to 20 mM, SDS to 0.5%, and proteinase K to 50 µg/ml final concentrations.
- b. Incubate at 65°C for 1 hr.
- c. Add an equal volume of phenol:chloroform, and mix by gentle inversion for 10 min.
- d. Centrifuge in a Beckman J2-21 centrifuge at 7,000 rpm for 10 min at room temperature (or for 5 min in a microcentrifuge for smaller volumes).
- e. Collect the supernatant, and repeat steps a–d until the interface is clean (usually one extraction is sufficient).
- f. Repeat steps a–d with chloroform only to remove any residual phenol (the small amount of aqueous material left at the interface can be collected separately in 1.5-ml microcentrifuge tubes, centrifuged, and combined with the larger fraction).
- g. Add $\frac{1}{10}$ volume of 3 M NaOAc and 2.5 volumes of 95% ethanol. Store at –20°C for at least 1 hr for DNA precipitation.
- h. Centrifuge the DNA at 12,000–15,000 rpm for 15 min at room temperature.
 - i. Pour off the supernatant, and wash the pellet with 70% ethanol.
 - j. Centrifuge for 5 min, and pour off the supernatant.

V. Library Protocols *continued*

- k. Spin down the residual ethanol. The last few drops can be removed with a micropipette.
- l. Let the pellet dry until the edges of the pellet begin to turn clear.

Notes:

- If the DNA is not completely dry, it will resuspend well in TE buffer. If the DNA is dry, add 1 ml of TE buffer, and allow the pellet to resuspend at 4°C overnight.
- The DNA should form a sharp band on an agarose gel with some contaminating RNA. RNase A can be added along with restriction enzymes, thus eliminating an additional digestion. A spectrophotometer reading should give characteristic curves, with a maximum absorption at 260 nm. The $A_{260/280}$ ratio should be >1.8.

(Alternative Protocol)

1. To 20 ml of phage liquid stock, add 4 µl of DNase I (5 mg/ml) and 75 µl of RNase A (10 mg/ml)
2. Add 20 ml of 20% PEG 8000/2M NaCl in 1X lambda dilution buffer. Incubate for 1 h at 0°C.
3. Centrifuge at 4°C for 20 min at 10,000 x g.
4. Remove supernatant by aspiration. Stand tube in an inverted position on a paper towel.
5. Add 0.5 ml of cold 1X lambda dilution buffer and resuspend by pipeting/scrapping off the pellet. Transfer to an Eppendorf tube.
6. Spin at 4°C for 2 min at 10,000 rpm. Transfer supernatant to a new tube.
7. Add 2.5 µl of 20% SDS, 5 µl of 0.5 M EDTA (pH 8), and 1.25 µl of 20mg/ml proteinase K (50 µg/ml final concentration).
8. DNA Extraction
 - a. Extract with 0.5 ml phenol. Vortex for 1 min, and centrifuge for 5 min at 12,500 rpm.
 - b. Extract with 0.5 ml phenol/chloroform. Vortex for 1 min, and centrifuge for 5 min at 12,500 rpm.
 - c. Extract with 0.5 ml chloroform. Vortex for 1 min, and centrifuge for 5 min at 12,500 rpm.
9. Precipitate
 - a. Add 1 ml of EtOH and mix by gentle inversion. DO NOT VORTEX. You should see a precipitate.
 - b. Spool out the DNA with a sealed capillary tube.
 - c. Wash with 70% EtOH.
10. Resuspend in 50 µl of TE. Use 2–5 µl to digest, and analyse on an agarose gel.

V. Library Protocols *continued*

G. Insert Excision

1. Restriction enzyme digestion

- Digest 5–10 μ l of purified DNA with the appropriate restriction enzyme at 37°C for 3–16 hr.
- Heat samples at 70°C for 10 min to inactivate the enzyme.
- Analyze the digest on a 0.8% agarose minigel.

Note: Depending on lysate titer, digestion of the entire isolate may be necessary in order to view the digest pattern. Precipitation of more than the usual volume of lysate can be performed (i.e., >10 ml) if a low titer is anticipated. If an abundance of undigested DNA is seen, the following steps are recommended.

- Reduce the DNA concentration.
- Increase the amount of restriction enzyme.
- Digest the DNA overnight.
- Combine the above suggestions.
- Purify the DNA using standard methods followed by a round of phage DNA purification on a CHROMA SPIN-400 Column.

2. Difficulty with insert excision

Inserts in both λ gt10 and λ gt11 libraries can be excised by EcoR I digestion. However, restriction enzymes are often sensitive to the purity of the DNA preparation. If the insert cannot be excised using EcoR I, the following options are available.

Note: Top agarose must be used during library plating, otherwise DNA purity is compromised and insert excision may prove difficult.

- Digest λ gt10 clones with Hind III/Ava I.

Digest λ gt11 clones with BsiW I.

Note: For additional restriction site information, please refer to Appendices B and C.

- 5'-STRETCH or 5'-STRETCH PLUS cDNA Libraries (5000 series)

The 5000 series libraries were constructed with a specially designed EcoR I adaptor that contains flanking Not I and Sal I sites. Therefore, the inserts can be excised with a Not I or Sal I digestion.

- NucleoBond®, NucleoSpin®, and NucleoTrap® Purification Products

Clontech offers a complete line of nucleic acid purification products that are useful for purifying library clones and inserts.

- Insert size determination by PCR[†]

Following identification of positive clones in λ gt10 or λ gt11, insert sizes can be determined by the use of Clontech's Insert Screening Amplimers. A listing of Clontech's cDNA Insert Screening Amplimers is presented in Section VII. Each Amplimer Set contains a complete protocol for use in PCR assays.

VI. Related Products

For DNA or oligonucleotide probe labeling	Cat. No.	Size
---	----------	------

- | | | |
|---------------------------|--------|--------|
| • Biotin-21-dUTP (0.5 mM) | 635701 | 100 µl |
| • Biotin-21-dUTP (10 mM) | 635702 | 100 µl |

Note: Clontech offers a wide range of products for nonisotopic labeling of probes. For additional information, please consult the current Clontech Catalog.

For library screening	Cat. No.	Size
-----------------------	----------	------

- | | | |
|--------------------------------------|--------|--------|
| • ExpressHyb™ Hybridization Solution | 636831 | 250 ml |
|--------------------------------------|--------|--------|

For phage DNA isolation and purification	Cat. No.	Size
--	----------	------

- | | | |
|---|------|--|
| • CHROMA SPIN™ Columns | many | |
| • NucleoBond®, NucleoSpin®, and NucleoTrap® Purification Products | many | |

VII. References

- Ausubel F.M. et al., eds. (1987–1995) *Current Protocols in Molecular Biology* (John Wiley & Sons, Inc.).
- Chenchik, A., Diatchenko, L., Chang, C. & Kuchibhatla, S. (January 1994) Great Lengths cDNA Synthesis Kit for high yields of full-length cDNA. *Clontechniques* **IX**(1):9–12.
- Dame, J. B., Williams, J. L., McCutchan, T. F., Weber, J. L., Wirtz, R. A., Hockmeyer, W. T., Maloy, W., L., Haynes, J. D., & Schneider, I. (1984) Structure of the gene encoding the immunodominant surface antigen on the sporozoite of the human malaria parasite *Plasmodium falciparum*. *Science* **225**:593–599.
- Gubler, U. & Hoffman, B. J. (1983) A simple and very efficient method for generating cDNA libraries. *Gene* **25**:263–269.
- Hagen, F. S., Gray, C. L. & Kuijper, J. L. (1988) Assaying the quality of cDNA libraries. *BioTechniques* **6**(4):340–345.
- Huynh, T. U., Young, R. A. & Davis, R. W. (1984) *DNACloning Techniques: A Practical Approach*, ed. Glover, D. (IRL Press, Oxford).
- Jacobs, K. A., Rudersdorf, R., Neill, S. D., Dougherty, J. P., Brown, E. L. & Fritsch, E. F. (1988) The thermal stability of oligonucleotide independent in tetraalkylammonium salt solutions: Application to identifying recombinant DNA clones. *Nucleic Acids Res.* **16**(10):4637–4650.
- Lathe, R. (1985) Synthetic oligonucleotide probes deduced from amino acid sequence data. Theoretical and practical considerations. *J. Mol. Biol.* **183**:1–12.
- Pirotta, V. (1979) *Nucleic Acids Res.* **6**:1495.
- Rosenberg, H. (1993) Bacterial chromosomal contamination of a eukaryotic cDNA library. *BioTechniques* **15**(2):200–201.
- Sambrook, J. & Russell, D.W. (2001) *Molecular Cloning: A Laboratory Manual*, (Cold Spring Harbor Laboratory Press, Cold Spring Harbor, NY).
- Suggs, S. V., Wallace, R. B., Hirose, T., Kawashima, E. H. & Itakura, K. (1981) Use of synthetic oligonucleotides as hybridization probes: Isolation of cloned cDNA sequences for human β 2-microglobulin. *Proc. Natl. Acad. Sci. USA* **78**:6613–6617.
- Wallace, R.B. et al. (1981) The use of synthetic oligonucleotides as hybridization probes II. Hybridization of oligonucleotides of mixed sequences to rabbit β -globin DNA. *Nucleic Acid Res.* **9**:879.
- Woo, S. L. C. (1979) A sensitive and rapid method for recombinant phage screening. *Methods Enzymol.* **68**:389.
- Wood, W. I., Gitschier, J., Lasky, L. A. & Lawn, R. M. (1985) Base composition-independent hybridization in tetramethylammonium chloride: A method for oligonucleotide screening of highly complex gene libraries. *Proc. Natl Acad. Sci. USA* **82**:1585–1588.
- Wood, W. B. (1966) Host specificity of DNA produced by *Escherichia coli*: Bacterial mutations affecting the restriction and modification of DNA. *J. Mol. Biol.* **16**:118.
- Young, R. A. & Davis, R. W. *Genetic Engineering*, eds. Setlow, J. & Hollaender (Plenum Press, NY).
- Young, R. A. & Davis, R. W. (1983) Efficient isolation of genes using antibody probes. *Proc. Natl. Acad. Sci. USA* **80**:1194–1198.
- Young, R. A. & Davis, R. W. (1983) Yeast RNA polymerase II genes: Isolation with antibody

Appendix A: Host Strain Information

BACTERIAL HOST STRAINS		
Strain	Reference	Genotype
C600	9–11	lacY1, leuB6, mcrA ⁻ , supE44, thi-1, thr-1, tonA21 For making lysates and for λgt10 propagation.
C600 hfl	10, 11	hflA150[chr::Tn10], lacY1, leuB6, mcrA ⁻ , supE44, thi-1, thr-1, tonA21 Note: Tn10 confers tet ^r For selection of λgt10 recombinants; the high frequency lysogeny mutation suppresses plaque formation by cl ⁺ bacteriophages while allowing plaque formation by recombinant cl ⁻ bacteriophages.
K802	12	gal K2, gal T22, hsdR2 (r _k ⁻ , m _k ⁺), lacY1, mcrA ⁻ , mcrB ⁻ , metB1, mrr ⁺ , supE44 For propagation of λ vectors and their recombinants.


Appendix A: Host Strain Information *continued*

BACTERIAL HOST STRAINS


Strain	Reference	Genotype
K803	12	gal K2, gal T22, hsdS3 (r_k^- , m_k^-), lacY1, mcrA ⁻ , mcrB ⁻ , metB1, mrr ⁺ , supE44 For propagation of λ vectors and their recombinants.
LE392	13, 14	gal K2, gal T22, hsdR514, lacY1, mcrA ⁻ , metB1, supE44, supF58, trpR55 For propagation of λ vectors and their recombinants; a derivative of ED8654.
Y1090r ⁻	11	araD139, hsdR (r_k^- , m_k^+), mcrA ⁻ , rpsL, supF, trpC22::Tn10, Δ lacU169, Δ lon, (pMC9) Note: rpsL confers str ^r , Tn10 confers tet ^r , and pMC9 confers tet ^r and amp ^r Strain of choice for immunological screening of expression libraries and propagation of λ gt11; expression of foreign proteins is controlled by high levels of lac repressor made by resident plasmid pMC9, which carries lacI ^q ; in lon ⁻ cells, β -galactosidase fusion proteins are often more stable than in wild-type cells; supF suppresses normally defective lysis of λ gt11 and λ gt22, leading to a high frequency of lytic plaques and easier antibody screening.

Appendix B: Vector Maps

FIGURE 4. λ gt10 λ gt11



Appendix B: Vector Maps *continued*



Appendix C: Sequence Information

1. λ gt10: Nucleic acid sequence around EcoR I Cloning site (15–17)

Note: EcoR I site is at position 251

```

1  v Hind III
AAGCTTCACA CCACGAACCA GCTCTAACCA TGCTAATCAA TGGATATTTC CCTTTGGGCT
TTCGAAGTGT GGTGCTTGGT CGAGATTGGT ACGATTAGTT ACCTATAAAG GGAAACCCGA

61                                     v BstB I
CAACGTGCCC AACAAATCTA ACATTCGAAT CAGAGGTGCC ATTGAGCAGC CAGTCAACAC
GTTGCACGGG TTGTTTAGAT TGTAAGCTTA GTCTCCACGG TAACTCGTCG GTCAGTTGTG

121
TTACGCCAAG AGCTGACGCA AGTTCCTGGTA AAAAGCGTGG TCGCTTAGTT TTACCGTTTT
AATGCGGTTT TCGACTGCGT TCAAGACCAT TTTTCGCACC AGCGAATCAA AATGGCAAAA

181
CGAGCTGCTC TATAGACTGC TGGGTAGTCC CCACCTTTTG AGCAAAGTTCA GCCTGGTTAA
GCTCGACGAG ATATCTGACG ACCCATCAGG GGTGGAAAAC TCGTTCAGT CGGACCAATT

241                                     v EcoR I
GTCCAAGCTG AATTC'TTTTG CTTTTTACCC TGGAAGAAAT ACTCATAAGC CACCTCTGTT
CAGGTTTCGAC TTAAGAAAAC GAAAAATGGG ACCTTCTTTA TGAGTATTCG GTGGAGACAA
                                     3' insert screening amplicer

301
ATTTACCCCC AATCTTCACA AGAAAACTG TATTTGACAA ACAAGATACA TTGTATGAAA
TAAATGGGGG TTAGAAGTGT TCTTTTTGAC ATAAACTGTT TGTCTATGT AACATACTTT

361                                     v Earl
ATACAAGAAA GTTTGTTGAT GGAGGCGATA TGCAAACCTT TCTGAACGC CTCAAGAAGA
TATGTCTCTT CAAACAACCT CCTCCGCTAT ACGTTTGAGA AAGACTTGCG GAGTTCTTCT

421
GGCGAATTGC GTTAAAAATG ACGCAAACCG AACTGGCAAC CAAAGCCGGT GTTAAACAGC
CCGCTTAACG CAATTTTAC TGCCTTTGGC TTGACCGTTG GTTTCGGCCA CAATTTGTGC

461
AATCAATTCA ACTGATTGAA
TTAGTTAAGT TGACTAACTT

```

IMM 434 cl region

5' insert screening amplicer: 5'-AGC AAG TTC AGCCTG GTT AAG T-3'

3' insert screening amplicer: 3'-GGG ACC TTC TTT ATG AGT ATT-5'

2. λ gt11: Nucleic acid sequence around Sac I–Kpn I region

Note: EcoR I site is at position 1070

```

1  v Sac I
GAGTCCTGC ACTGGATGGT GGCCTGGAT GGTAAGCCGC TGGCAAGCGG TGAAGTGCCT
CTCGAGGACG TGACCTACCA CCGCGACCTA CCATTTCGGC ACCGTTTCGC ACTTCACGGA

```

Appendix C: Sequence Information *continued*

61	CTGGATGTCG	CTCCACAAGG	TAAACAGTTG	ATTGAACTGC	CTGAACCTACC	GCAGCCGGAG
	GACCTACAGC	GAGGTGTTC	ATTTGTCAAC	TAACCTTGACG	GACTTGATGG	CGTCGGCCTC
121	AGCGCCGGGC	AACTCTGGCT	CACAGTACGC	GTAGTGCAAC	CGAACGCGAC	CGCATGGTCA
	TCGCGGCCCG	TTGAGACCGA	GTGTCATGCG	CATCACGTTG	GCTTGCGCTG	GCGTACCAGT
181	GAAGCCGGGC	ACATCAGCGC	CTGGCAGCAG	TGGCGTCTGG	CGGAAAACCT	CAGTGTGACG
	CTTCGGCCCCG	TGTAGTCGCG	GACCGTCGTC	ACCGCAGACC	GCCTTTTGA	GTCACACTGC
241	CTCCCCGCGG	CGTCCCACGC	CATCCCGCAT	CTGACCACCA	GCGAAATGGA	TTTTTGTCATC
	GAGGGGCGGC	GCAGGGTGCG	GTAGGGCGTA	GACTGGTGGT	CGCTTTACCT	AAAAACGTAG
301	GAGCTGGGTA	ATAAGCGTTG	GCAATTTAAC	CGCCAGTCAG	GCTTCTTTTC	ACAGATGTGG
	CTCGACCCAT	TATTCGCAAC	CGTTAAATTG	GCGGTCAGTC	CGAAAGAAAG	TGTCTACACC
361	ATTGGCGATA	AAAAACAAC	GCTGACGCCG	CTGCGCGATC	AGTTCACCCG	TGCACCCGTC
	TAACCGCTAT	TTTTTGTGA	CGACTGCGGC	GACGCGCTAG	TCAAGTGGGC	ACGTGGCGAC
421	GATAACGACA	TTGGCGTAAG	TGAAGCGACC	CGCATTGACC	CTAACGCCTG	GGTCGAACGC
	CTATTGCTGT	AACCGCATTC	ACTTCGCTGG	GCGTAACCTG	GATTGCGGAC	CCAGCTTGCG
481	TGGAAGGCGG	CGGGCCATTA	CCAGGCCGAA	GCAGCGTTGT	TGCAGTGCAC	GGCAGATACA
	ACCTTCCGCC	GCCCGGTAAT	GGTCCGGCTT	CGTCGCAACA	ACGTACAGTG	CCGTCTATGT
541	CTTGCTGATG	CGGTGCTGAT	TACGACCGCT	CACGCGTGCG	AGCATCAGGG	GAAAACCTTA
	GAACGACTAC	GCCACGACTA	ATGCTGGCGA	GTGCGCACCG	TCGTAGTCCC	CTTTTGGAAT
601	TTTATCAGCC	GGAAAACCTA	CCGATTGAT	GGTAGTGCTC	AAATGGCGAT	TACCGTTGAT
	AAATAGTCGG	CCTTTTGAT	GGCCTAACTA	CCATCACCAG	TTTACCGCTA	ATGGCAACTA
661	GTTGAAGTGG	CGAGCGATAC	ACCGCATCCG	GCGCGGATTG	GCCTGAACCTG	CCAGCTGGCG
	CAACTTCACC	GCTCGCTATG	TGGCGTAGGC	GCGCCCTAAC	CGGACTTGAC	GGTCGACCGC
721	CAGGTAGCAG	AGCGGGTAAA	CTGGCTCGGA	TTAGGGCCGC	AAGAAAACCTA	TCCCGACCGC
	GTCCATCGTC	TCGCCCATTT	GACCGAGCCT	AATCCCGCGG	TTCTTTTGTAT	AGGGCTGGCG
781	CTTACTGCCG	CCTGTTTTGA	CCGCTGGGAT	CTGCCATTGT	CAGACATGTA	TACCCCGTAC
	GAATGACGGC	GGACAAAAC	GGCGACCTTA	GACGGTAACA	GTCTGTACAT	ATGGGGCATG

Appendix C: Sequence Information *continued*

841	GTCTTCCCGA	GCGAAAACGG	TCTGCGCTGC	GGGACGCGCG	AATTGAATTA	TGGCCACAC
	CAGAAGGGCT	CGCTTTTGCC	AGACGCGACG	CCCTGCGCGC	TTAACTTAAT	ACCGGGTGTG
901	ν BstX I	ν RleA			ν AlwN I	
	CAGTGGCGCG	GCGACTTCCA	GTTCAACATC	AGCCGCTACA	GTCAACAGCA	ACTGATGGAA
	GTCACCGCGC	CGCTGAAGGT	CAAGTTGTAG	TCGGCGATGT	CAGTTGTCGT	TGACTACCTT
961	ν Bsg I					
	ACCAGCCATC	GCCATCTGCT	GCACGCGGAA	GAAGGCACAT	GGCTGAATAT	CGACGGTTTC
	TGGTCGGTAG	CGGTAGACGA	CGTGCGCCTT	CTTCGGTGTA	CCGACTTATA	GCTGCCAAAG
1021	ν Nde I ν Xcm I		5' insert screening amplimer		Gsu I ν EcoR I	Esp I ν
	CATATGGGGA	TTGGTGGCGA	CGACTCCTGG	AGCCCGTCAG	TATCGGCGGA	ATTCCAGCTG
	GTATACCCCT	AACCACCGCT	GCTGAGGACC	TCGGGCAGTC	ATAGCCGCCT	TAAGGTCGAC
1081						
	AGCGCCGGTC	GCTACCATTA	CCAGTTGGTC	TGGTGTCAAA	AATAATAATA	ACCGGGCAGG
	TCGCGGCCAG	CGATGGTAAT	GGTCAACCAG	ACCACAGTTT	TTATTATTAT	TGGCCCGTCC
	3' insert screening amplimer					
1141						
	CCATGTCTGC	CCGTATTTTCG	CGTAAGGAAA	TCCATTGTAC	TGCCGGACCA	CCGACTGTGA
	GGTACAGACG	GGCATAAAGC	GCATTCTCTT	AGGTAACATG	ACGGCTGGT	GGCTGACACT
1201		Xcm I				
	Bgl I ν ν	ν Nco I	ν BsiW I	ν BspM I	Age I ν	ν Rsr II
	GCCACTCCGG	CCATGGCGTA	CGCACTGACC	TGCTTACTGA	TTTGTAAGAC	CGGTCCGGCC
	CGGTGAGGCC	GGTACCGCAT	GCGTGACTGG	ACGAATGACT	AAACATTTTG	GCCAGGCCGG
1261						
	ATCACGCTCA	CATAACGTCC	ACGCAGGCTC	TCATAGTGAA	ACGTATCCTC	CCCGGTCATC
	TAGTGCGAGT	GTATTGCAGG	TGCGTCCGAG	AGTATCACTT	TGCATAGGAG	GGGCCAGTAG
1321						
	ACTGTGCTGC	TCTTTTTTCGA	CGCGGCAGAC	CCCAGGGAAG	CCATCACCCC	CACACTGTCC
	TGACACGACG	AGAAAAAGCT	GCGCCGCTTG	GGGTCCCTTC	GGTAGTGGGG	GTGTGACAGG
1381		ν RleA I			ν Age I	
	GTCAGTCTAT	AACGGTACTT	CACGTTAATC	CCTTTCAGAT	GACTCACACC	GGTATCCCCG
	CAGTCGAGTA	TTGCCATGAA	GTGCAATTAG	GGAAAGTCTA	CTGAGTGTGG	CCATAGGGGC
1441						
	CCCACAACG	ACGGCAATGT	ACCCGGTTTC	ACTTGAAAAT	AGCCACCGT	AAACGTACCA
	GGGCTGTTGC	TGCCGTTACA	TGGGCCAAAG	TGAACTTTTA	TCGGGTGGCA	TTTGCATGGT
1501						
	TGTCCACCTT	CCGCACGGGC	CGGAGTGACT	GTCACCGCAA	GTGCGGCAAA	GACAGCAACG
	ACAGGTGGAA	GGCGTGCCCG	GCCTCACTGA	CAGTGCGGTT	CACGCCGTTT	CTGTGCTTGC
1561						
	GCAATACACA	CATTACGCAT	CGTTCACCTC	TCACTGTTTT	ATAATAAAAC	GCCCGTCCCC
	CGTTATGTGT	GTAATGCGTA	GCAAGTGAG	AGTGACAAAA	TATTATTTTG	CGGGCAAGGG

Appendix C: Sequence Information *continued*

There are no restriction sites in the previous sequence for the following enzymes:

Aat II	Afl II	ApaB I	Apa I
Asc I	Ase I	Avr II	BamH I
Bbe I	Bcg I	Bcl I	Bgl II
Bpu 10 I	Bsa I	Bsi I	Bsp120 I
BspE I	BspH I	BssH II	BstB I
BstE II	Bsu361	Cla I	Dra III
Eag I	Ear I	Eco571	EcoN I
Esp3 I	Fse I	HgiEI I	Hind III
Kas I	Mun I	Nar I	NgoM I
Nhe I	Not I	Nsi I	Pac I
PpuM I	Pst I	Pvu I	Sac II
Sal I	Sap I	Sce I	Sfi I
SgrA I	Spe I	Sph I	Sse I
TthIII I	Xba I	Xho I	Xma I

3. λ gt11: Amino acid sequence around Sac I–Kpn I region

1
GAG CTC CTG CAC TGG ATG GTG GCG CTG GAT GGT AAG CCG CTG GCA AGC GGT GAA
CTC GAG GAC GTG ACC TAC CAC CGC GAC CTA CCA TTC GGC GAC CGT TCG CCA CTT
Glu Leu Leu His Trp Met Val Ala Leu Asp Gly Lys Pro Leu Ala Ser Gly Glu

55
GTG CCT CTG GAT GTC GCT CCA CAA GGT AAA CAG TTG ATT GAA CTG CCT GAA CTA
CAC GGA GAC CTA CAG CGA GGT GTT CCA TTT GTC AAC TAA CTT GAC GGA CTT GAT
Val Pro Leu Asp Val Ala Pro Gln Gly Lys Gln Leu Ile Glu Leu Pro Glu Leu

110
CCG CAG CCG GAG AGC GCC GGG CAA CTC TGG CTC ACA GTA CGC GTA GTG CAA CCG
GGC GTC GGC CTC TCG CGG CCC GTT GAG ACC GAG TGT CAT GCG CAT CAC GTT GGC
Pro Gln Pro Glu Ser Ala Gly Gln Leu Trp Leu Thr Val Arg Val Val Gln Pro

164
AAC GCG ACC GCA TGG TCA GAA GCC GGG CAC ATC AGC GCC TGG CAG CAG TGG CGT
TTG CGC TGG CGT ACC AGT CTT CGG CCC GTG TAG TCG CGG ACC GTC GTC ACC GCA
Asn Ala Thr Ala Trp Ser Glu Ala Gly His Ile Ser Ala Trp Gln Gln Trp Arg

218
CTG GCG GAA AAC CTC AGT GTG ACG CTC CCC GCC GCG TCC CAC GCC ATC CCG CAT
GAC CGC CTT TTG GAG TCA CAC TGC GAG GGG CGG CGC AGG GTG CGG TAG GGC GTA
Leu Ala Glu Asn Leu Ser Val Thr Leu Pro Ala Ala Ser His Ala Ile Pro His

272
CTG ACC ACC AGC GAA ATG GAT TTT TGC ATC GAG CTG GGT AAT AAG CGT TGG CAA
GAC TGG TGG TCG CTT TAC CTA AAA ACG TAG CTC GAC CCA TTA TTC GCA ACC GTT
Leu Thr Thr Ser Glu Met Asp Phe Cys Ile Glu Leu Gly Asn Lys Arg Trp Gln

Appendix C: Sequence Information *continued*

326

TTT AAC CGC CAG TCA GGC TTT CTT TCA CAG ATG TGG ATT GGC GAT AAA AAA CAA
AAA TTG GCG GTC AGT CCG AAA GAA AGT GTC TAC ACC TAA CCG CTA TTT TTT GTT
Phe Asn Arg Gln Ser Gly Phe Leu Ser Gln Met Trp Ile Gly Asp Lys Lys Gln

380

CTG CTG ACG CCG CTG CGC GAT CAG TTC ACC CGT GCA CCG CTG GAT AAC GAC ATT
GAC GAC TGC GGC GAC GCG CTA GTC AAG TGG GCA CGT GGC GAC CTA TTG CTG TAA
Leu Leu Thr Pro Leu Arg Asp Gln Phe Thr Arg Ala Pro Leu Asp Asn Asp Ile

434

GGC GTA AGT GAA GCG ACC CGC ATT GAC CCT AAC GCC TGG GTC GAA CGC TGG AAG
CCG CAT TCA CTT CGC TGG GCG TAA CTG GGA TTG CGG ACC CAG CTT GCG ACC TTC
Gly Val Ser Glu Ala Thr Arg Ile Asp Pro Asn Ala Trp Val Glu Arg Trp Lys

488

GCG GCG GGC CAT TAC CAG GCC GAA GCA GCG TTG TTG CAG TGC ACG GCA GAT ACA
CGC CGC CCG GTA ATG GTC CGG CTT CGT CGC AAC AAC GTC ACG TGC CGT CTA TGT
Ala Ala Gly His Tyr Gln Ala Glu Ala Ala Leu Leu Gln Cys Thr Ala Asp Thr

542

CTT GCT GAT GCG GTG CTG ATT ACG ACC GCT CAC GCG TGG CAG CAT CAG GGG AAA
GAA CGA CTA CGC CAC GAC TAA TGC TGG CGA GTG CGC ACC GTC GTA GTC CCC TTT
Leu Ala Asp Ala Val Leu Ile Thr Thr Ala His Ala Trp Gln His Gln Gly Lys

596

ACC TTA TTT ATC AGC CGG AAA ACC TAC CGG ATT GAT GGT AGT GGT CAA ATG GCG
TGG AAT AAA TAG TCG GCC TTT TGG ATG GCC TAA CTA CCA TCA CCA GTT TAC CCG
Thr Leu Phe Ile Ser Arg Lys Thr Tyr Arg Ile Asp Gly Ser Gly Gln Met Ala

650

ATT ACC GTT GAT GTT GAA GTG GCG AGC GAT ACA CCG CAT CCG GCG CGG ATT GGC
TAA TGG CAA CTA CAA CTT CAC CGC TCG CTA TGT GGC GTA GGC CGC GCC TAA CCG
Ile Thr Val Asp Val Glu Val Ala Ser Asp Thr Pro His Pro Ala Arg Ile Gly

704

CTG AAC TGC CAG CTG GCG CAG GTA GCA GAG CGG GTA AAC TGG CTC GGA TTA GGG
GAC TTG ACG GTC GAC CGC GTC CAT CGT CTC GCC CAT TTG ACC GAG CCT AAT CCC
Leu Asn Cys Gln Leu Ala Gln Val Ala Glu Arg Val Asn Trp Leu Gly Leu Gly

758

CCG CAA GAA AAC TAT CCC GAC CGC CTT ACT GCC GCC TGT TTT GAC CGC TGG GAT
GGC GTT CTT TTG ATA GGG CTG GCG GAA TGA CGG CGG ACA AAA CTG GCG ACC CTA
Pro Gln Glu Asn Tyr Pro Asp Arg Leu Thr Ala Ala Cys Phe Asp Arg Trp Asp

812

CTG CCA TTG TCA GAC ATG TAT ACC CCG TAC GTC TTC CCG AGC GAA AAC GGT CTG
GAC GGT AAC AGT CTG TAC ATA TGG GGC ATG CAG AAG GGC TCG CTT TTG CCA GAC
Leu Pro Leu Ser Asp Met Tyr Thr Pro Tyr Val Phe Pro Ser Glu Asn Gly Leu

Appendix C: Sequence Information *continued*

866

CGC TGC GGG ACG CGC GAA TTG AAT TAT GGC CCA CAC CAG TGG CGC GGC GAC TTC
GCG ACG CCC TGC GCG CTT AAC TTA ATA CCG GGT GTG GTC ACC GCG CCG CTG AAG
Arg Cys Gly Thr Arg Glu Leu Asn Tyr Gly Pro His Gln Trp Arg Gly Asp Phe

920

CAG TTC AAC ATC AGC CGC TAC AGT CAA CAG CAA CTG ATG GAA ACC AGC CAT CGC
GTC AAG TTG TAG TCG GCG ATG TCA GTT GTC GTT GAC TAC CTT TGG TCG GTA GCG
Gln Phe Asn Ile Ser Arg Tyr Ser Gln Gln Gln Leu Met Glu Thr Ser His Arg

974

CAT CTG CTG CAC GCG GAA GAA GGC ACA TGG CTG AAT ATC GAC GGT TTC CAT ATG
GTA GAC GAC GTG CGC CTT CTT CCG TGT ACC GAC TTA TAG CTG CCA AAG GTA TAC
His Leu Leu His Ala Glu Glu Gly Thr Trp Leu Asn Ile Asp Gly Phe His Met

1028 EcoR I

GGG ATT GGT GGC GAC GAC TCC TGG AGC CCG TCA GTA TCG GCG **GAA TTC** CAG CTG
CCC TAA CCA CCG CTG CTG AGG ACC TCG GGC AGT CAT AGC CGC CTT AAG GTC GAC
Gly Ile Gly Gly Asp Asp Ser Trp Ser Pro Ser Val Ser Ala Glu Phe Gln Leu

1082

AGC GCC GGT CGC TAC CAT TAC CAG TTG GTC TGG TGT CAA AAA TAA TAA TAA CCG
TCG CGG CCA GCG ATG GTA ATG GTC AAC CAG ACC ACA GTT TTT ATT ATT ATT GGC
Ser Ala Gly Arg Tyr His Tyr Gln Leu Val Trp Cys Gln Lys End End End Pro

1136

GGC AGG CCA TGT CTG CCC GTA TTT CGC GTA AGG AAA TCC ATT GTA CTG CCG GAC
CCG TCC GGT ACA GAC GGG CAT AAA GCG CAT TCC TTT AGG TAA CAT GAC GGC CTG
Gly Arg Pro Cys Leu Pro Val Phe Arg Val Arg Lys Ser Ile Val Leu Pro Asp

1190

CAC CGA CTG TGA GCC ACT CCG GCC ATG GCG TAC GCA CTG ACC TGC TTA CTG ATT
GTG GCT GAC ACT CGG TGA GGC CGG TAC CGC ATG CGT GAC TGG ACG AAT GAC TAA
His Arg Leu End Ala Thr Pro Ala Met Ala Tyr Ala Leu Thr Cys Leu Leu Ile

1244

TGT AAA ACC GGT CCG GCC ATC ACG CTC ACA TAA CGT CCA CGC AGG CTC TCA TAG
ACA TTT TGG CCA GGC CGG TAG TGC GAG TGT ATT GCA GGT GCG TCC GAG AGT ATC
Cys Lys Thr Gly Pro Ala Ile Thr Leu Thr End Arg Pro Arg Arg Leu Ser End

1298

TGA AAC GTA TCC TCC CCG GTC ATC ACT GTG CTG CTC TTT TTC GAC GCG GCG AAC
ACT TTG CAT AGG AGG GGC CAG TAG TGA CAC GAC GAG AAA AAG CTG CGC CGC TTG
End Asn Val Ser Ser Pro Val Ile Thr Val Leu Leu Phe Phe Asp Ala Ala Asn

1352

CCC AGG GAA GCC ATC ACC CCC ACA CTG TCC GTC AGC TCA TAA CGG TAC TTC ACG
GGG TCC CTT CGG TAG TGG GGG TGT GAC AGG CAG TCG AGT ATT GCC ATG AAG TGC
Pro Arg Glu Ala Ile Thr Pro Thr Leu Ser Val Ser Ser End Arg Tyr Phe Thr

Appendix C: Sequence Information *continued*

1406

TTA ATC CCT TTC AGA TGA CTC ACA CCG GTA TCC CCG CCC GAC AAC GAC GGC AAT
AAT TAG GGA AAG TCT ACT GAG TGT GGC CAT AGG GGC GGG CTG TTG CTG CCG TTA
Leu Ile Pro Phe Arg End Leu Thr Pro Val Ser Pro Pro Asp Asn Asp Gly Asn

1460

GTA CCC GGT TTC ACT TGA AAA TAG CCC ACC GTA AAC GTA CCA TGT CCA CCT TCC
CAT GGG CCA AAG TGA ACT TTT ATC GGG TGG CAT TTG CAT GGT ACA GGT GGA AGG
Val Pro Gly Phe Thr End Lys End Pro Thr Val Asn Val Pro Cys Pro Pro Ser

1514

GCA CGG GCC GGA GTG ACT GTC ACC GCA AGT GCG GCA AAG ACA GCA ACG GCA ATA
CGT GCC CGG CCT CAC TGA CAG TGG CGT TCA CGC CGT TTC TGT CGT TGC CGT TAT
Ala Arg Ala Gly Val Thr Val Thr Ala Ser Ala Ala Lys Thr Ala Thr Ala Ile

1568

CAC ACA TTA CGC ATC GTT CAC CTC TCA CTG TTT TAT AAT AAA ACG CCC GTT CCC
GTG TGT AAT GCG TAG CAA GTG GAG AGT GAC AAA ATA TTA TTT TGC GGG CAA GGG
His Thr Leu Arg Ile Val His Leu Ser Leu Phe Tyr Asn Lys Thr Pro Val Pro

1622

GGA CGA ACC TCT GTA ACA CAC TCA GAC CAC GCT GAT GCC CAG CGC CTG TTT CTT
CCT GCT TGG AGA CAT TGT GTG AGT CTG GTG CGA CTA CGG GTC GCG GAC AAA GAA
Gly Arg Thr Ser Val Thr His Ser Asp His Ala Asp Ala Gln Arg Leu Phe Leu

1676

AAT CAC CAT AAC CTG CAC ATC GCT GGC AAA CGT ATA CGG CGG AAT ATC TGC CGA
TTA GTG GTA TTG GAC GTG TAG CGA CCG TTT GCA TAT GCC GCC TTA TAG ACG GCT
Asn His His Asn Leu His Ile Ala Gly Lys Arg Ile Arg Arg Asn Ile Cys Arg

1730

ATG CCG TGT GGA CGT AAG CGT GAA CGT CAG GAT CAC GTT TCC CCG ACC CGC TGG
TAC GGC ACA CCT GCA TTC GCA CTT GCA GTC CTA GTG CAA AGG GGC TGG GCG ACC
Met Pro Cys Gly Arg Lys Arg Glu Arg Gln Asp His Val Ser Pro Thr Arg Trp

1784

CAT GTC AAC AAT ACG GGA GAA CAC CTG TAC CGC CTC GTT CGC GCG GCC ATC ATA
GTA CAG TTG TTA TGC CCT CTT GTG GAC ATG GCG GAG CAA GCG GCG CGG TAG TAT
His Val Asn Asn Thr Gly Glu His Leu Tyr Arg Leu Val Arg Arg Ala Ile Ile

1838

AAT CAC CGC ACC GTT CAT CAG TAC TTT CAG ATA ACA CAT CGA ATA CGT TGT CCT
TTA GTG GCG TGG CAA GTA GTC ATG AAA GTC TAT TGT GTA GCT TAT GCA ACA GGA
Asn His Arg Thr Val His Gln Tyr Phe Gln Ile Thr His Arg Ile Arg Cys Pro

1892

GCC GCT GAC AGT ACG CTT ACT TCC GCG AAA CGT CAG CGG AAG CAC CAC TAT CTG
CGG CGA CTG TCA TGC GAA TGA AGG CGC TTT GCA GTC GCC TTC GTG GTG ATA GAC
Ala Ala Asp Ser Thr Leu Thr Ser Ala Lys Arg Gln Arg Lys His His Tyr Leu

Appendix C: Sequence Information *continued*

1946

GCG ATC AAA AGG ATG GTC ATC GGT CAC GGT GAC AGT ACG GGT ACC
CGC TAG TTT TCC TAC CAG TAG CCA GTG CCA CTG TCA TGC CCA TGG
Ala Ile Lys Arg Met Val Ile Gly His Gly Asp Ser Thr Gly Thr

Appendix D: Recombination Frequency Determination

The following protocol is used to determine the recombination frequency for libraries cloned in λ expression vectors (e.g., λ gt11).

Prepare recommended titering dilutions of phage lysate (library) as follows:

- 1. Pipet 2 μ l of the library lysate into 1 ml of 1X lambda dilution buffer (Dilution 1 = 1:500).
- 2. Transfer 2 μ l of Dilution 1 into a second tube containing 1 ml of 1X lambda dilution buffer (Dilution 2 = 1:250,000).
- 3. Prepare 4 tubes as described in Table IV using the bacterial overnight culture obtained in Step VI.B.1.

TABLE IV: LIBRARY PLATING DILUTIONS

Tube	1X Lambda Dilution Buffer	Bacterial Overnight Culture	Phage Dilution 2
1	100 μ l	200 μ l	2 μ l
2	100 μ l	200 μ l	5 μ l
3	100 μ l	200 μ l	10 μ l
4 (Control)	100 μ l	200 μ l	0 μ l

- 4. Incubate tubes in a 37°C water bath for 15 min.
- 5. Add 30 μ l of 0.1 M X-Gal and 30 μ l of 0.1 M IPTG to each of the 4 tubes.
- 6. Add 3 ml of melted LB soft top agar + MgSO₄ (45°C) to each of the 4 tubes. Mix well.
- 7. Pour the contents from each tube onto 4 separate 90-mm LB agar + MgSO₄ plates. Swirl the plates quickly after pouring to allow for even spreading of the agar.

Note: If the plates were prepared fresh and have been sitting at room temperature, prewarm plates at 37°C for 10–15 min. If plates have been stored at 4°C, prewarm plates at 37°C for 1 hr prior to use.

- 8. Cool plates at room temperature for 10 min to allow the inoculum to soak into the agar.

Appendix D: Recombination Frequency *continued*

9. Incubate plates (inverted position) at 37°C for at least 6–7 hr.
10. Count the number of blue plaques and the number of white plaques.

$$\text{Recombination frequency} = \frac{\text{number of blue plaques}}{\text{number of total plaques}} \times 100\%$$

Appendix E: rpm Calculation

$$\begin{aligned} \text{rpm} &= 1000 \times \sqrt{\text{RCF}/(1.12r)} \\ r &= \text{radius of rotor (in mm)} \\ \text{RCF} &= \text{relative centrifugal field} \\ g \text{ Force} &= \text{RCF} \times g \end{aligned}$$

Examples:

- a. Swinging bucket rotor = Beckman Model TJ-6

r	=	125 mm
RCF	=	700
therefore, rpm	=	2236
- b. Fixed-angle rotor = Eppendorf Model 5415C (microcentrifuge)

r	=	50 mm
RCF	=	700
therefore, rpm	=	3536

General Notice to Purchaser

Clontech products are to be used for research purposes only. They may not to be used for any other purpose, including, but not limited to, use in drugs, in vitro diagnostic purposes, therapeutics, or in humans. Clontech products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without written approval of Clontech Laboratories, Inc.

UltraClear™ is a trademark of Beckman Instruments.

† The Polymerase Chain Reaction (PCR) process is covered by U.S. Patents owned by Hoffmann-LaRoche Inc., and pending and issued patents owned by Hoffman-LaRoche Inc. in non-U.S. countries.

NucleoBond®, NucleoTrap®, and NucleoSpin® are registered trademarks of Macherey-Nagel GmbH & Co., KG.

Clontech, the Clontech logo and all other trademarks are the property of Clontech Laboratories, Inc. Clontech is a Takara Bio Company. ©2006