

Oligonucleotides

For Life Science Research

- > Custom **oligonucleotides**
- > **NGS** oligonucleotides
- > **RNAi** oligonucleotides
- > **Real-Time qPCR** probes
- > **Highly complex** oligonucleotides
- > **Catalogue** oligonucleotides
- > **Synthesis** reagents





A RELIABLE EXPERIENCE

Since 1985 Eurogentec has provided high-quality reagents and custom-synthesised oligonucleotides to scientists around the globe.


CUSTOM PRODUCTS

- All types of oligonucleotides, from 2 to 225 bases
- All chemistries: DNA, RNA, LNA®, 2'O-Me, 2'O-MOE, PNA,...
- More than 300 modifications
- All synthesis scales, from µg to grams
- Wide range of Real-Time qPCR Probes
- RNAi oligonucleotides
- Custom fill & finish

TRUSTED QUALITY

- Optimised chemistry
- Stringent quality controls
 - ISO 9001 certified quality system
 - ISO 13485 certified for IVD oligonucleotides


Research, Track™ and GMP grade oligonucleotides



From lab scale to large synthesis scale



More than 300 modifications with a wide range of dyes and quenchers



ISO 9001
ISO 13485 certified quality systems



Oli&GO™ oligonucleotide e-commerce platform



Life Science Oligonucleotides

Custom

Catalogue

Custom oligos

NGS

RNAi

qPCR probes

Unique™ oligos

Catalogue oligos

Synthesis reagents

- p5 *A large range for any applications*
- p9 **Custom oligos**
- p10 *NGS oligos*
- p11 *RNAi oligos*
- p13 *qPCR probes*
- p16 *Spectral properties of fluorophores & quenchers*
- p18 *Unique™ oligos*
- p19 **Catalogue oligos**
- p21 *Synthesis reagents*
- p22 *Additional services*

Annexes

- p23 *Synthesis scale vs guaranteed yield*
- p24 *Shipping*
- p24 *Documentation*
- p25 *How to store your oligo*
- p25 *How to reconstitute your oligo*
- p25 *How to quantify your oligo*
- p26 *IUB code*
- p27 *Dyes compatibility table*
- p28 *How to order*
- p28 *How to pay*
- p29 *How to reduce my shipping fees*
- p29 *Related products*
- p30 *License statements*
- p31 *Trademarks and labels*

A LARGE RANGE FOR ANY APPLICATION

Whatever your application, even for those that are most demanding (NMR, X-ray crystallography, *in vivo* animal studies...), Eurogentec can provide the highest quality oligonucleotides to meet (and exceed) your expectations!

CUSTOM OLIGOS p.9

PCR | FISH | Pyrosequencing
| Cloning | NMR | X-Ray
crystallography | Mutagenesis |
SNP Analysis

NGS OLIGOS p.10

Next-generation sequencing

RNAi OLIGOS p.11

Gene silencing | Antisense studies

qPCR PROBES p.13

Real-time qPCR | Patient
management | Diagnostic assays

UNIQUE™ OLIGOS p.18

For highly complex oligos or if
you don't find what you need,
please contact us at unique@eurogentec.com

CATALOGUE OLIGOS p.19


Cloning | Sequencing | PNA FISH |
qPCR Calibration




CHOOSE THE RECOMMENDED PURIFICATION

The aim of any purification method is to remove the by-products resulting from the removal of the protecting groups and other synthesis by-products. To know the best purification according to each modification, consult the price list available on www.eurogentec.com. If you are not sure which purification suits your application, then please specify "Recommended Single Purification" (additional fee) and we will choose the best purification for you.

PURIFICATION VS APPLICATIONS



PURIFICATION VS OLIGO LENGTH



*These values are purely indicative and only valid for an unmodified oligonucleotide of 20 bases. In addition, according to your oligonucleotides (sequences, modifications...), the purity level can be analysed by various methods (analytical HPLC, CGE...).

EUROGENTEC MANUFACTURES HIGHLY PURIFIED OLIGONUCLEOTIDES UP TO 95%.

SYNTHESIS SCALE VS GUARANTEED YIELD

Please refer to the minimum guaranteed yield table page 23 to select the right synthesis scale or contact us at: oligocentre@eurogentec.com

Custom oligos

Eurogentec proposes a large choice of chemistries, modifications, specifications and purifications. More than 300 modifications and several purity levels are available. ■

SPECIFICATIONS

- Length:** From 5 to 139 bases
- Synthesis scale:** 10 nmol • 40 nmol • 200 nmol • 1000 nmol • 2.5 μmol • 5 μmol • 10 μmol*
- Backbone:** DNA, RNA, LNA®, 2'O-Me RNA, 2'O-MOE RNA, PNA and all linkages
- Modifications:** More than 300 modifications! (see p. 7)
- Purifications:** SePOP desalting, RP-Cartridge•Gold™, HPLC, PAGE, Dual HPLC, UltraPureGold™
- Quality Control:** MALDI-TOF MS
- Format:** Dried (except for unmodified SePOP desalted oligonucleotides from 15 to 39 DNA bases: 100 μM H₂O by default)
- Packaging:** 2 mL tube, 96-well or 384-well plates
- Documentation:** Technical data sheet
- Shipping:** At room temperature


DID YOU KNOW?

The combination of our vast expertise in oligonucleotide to the well reputed know how of AnaSpec in complex peptide synthesis, allows us to offer you high quality **peptide-oligo conjugates**.

SPECIFIC NEED?

Need a high production process traceability? Discover our Track oligos on www.eurogentec.com/track-oligos.html

>Note
* Larger synthesis scales are available on request.



NGS oligos

Next-Generation Sequencing (NGS) is a high-throughput technology allowing the **massive sequencing** of nucleic acids following a DNA library preparation.

After DNA fragmentation, adapters (including indexes for multiplexing) are fused to the fragments. NGS adapters require both a **high level of purity** (no n-x side products) and the **absence of cross-contamination** (confusing index sequences).

Thanks to our long history as an oligo provider Eurogentec has developed a dedicated manufacturing process for the production of high quality NGS oligos. ■

SPECIFICATIONS

Quality: Low cross-contamination (<0,1%)

Length: from 20 to 85 bases

Quantity: 10 nmol minimum delivered*

Purification: HPLC or cartridge

QC: 100% QC checked by Maldi-TOF MS

5' Modifications: 5' Phosphate / 5' Biotin-TEG

Bases Option: Phosphorothioate bond

Wobble Bases: Available at no additional cost

Format: dried in tubes

Free shipping

RNAi oligos

WHAT IS RNAi INTERFERENCE?

RNA interference is a mechanism of gene silencing at the mRNA level. This phenomenon is triggered by small interfering (si)RNAs and micro (mi)RNAs.

siRNAs and miRNAs regulate gene expression. They can activate the degradation of the targeted mRNA or prevent its translation. ■

SPECIFICATIONS

Length: From 21 to 27 bases

Synthesis scale: 10 nmol · 40 nmol · 200 nmol · 1000 nmol*

Backbone: RNA, LNA®, 2'O-Me RNA, 2'O-MOE RNA and all linkages

Modifications: 5':Phosphate, 6-FAM, Cy®3, Cy5®, TET, HEX,...
3': DABCYL, TEG-Cholesteryl, TAMRA...

Purifications: SePOP Desalting or IEX-RP/HPLC

Quality Control: MALDI-TOF MS


Format: Dried

Packaging: 2 mL tube

Documentation: Technical data sheet

siRNA Design: Free and guaranteed

Shipping: At room temperature



DID YOU KNOW?

Eurogentec can synthesise highly modified and very stable **RNA oligos**.
Contact us at:
unique@eurogentec.com

Design assistance

info@eurogentec.com

Eurogentec all along your NGS process

EUROGENTEC
NGS
Oligonucleotides

Library Preparation

Amplification

Sequencing

Data Analysis
Alignment

>Note

*Larger amounts are available on request

>Note

- The antisense strand must either have a free 5'-OH (by default) or 5'-phosphate terminus.

- Certain modifications can sometimes be useful to increase stability or cellular uptake e.g. Modifying siRNA with cholesterol is used to facilitate tissue / cellular uptake.

- Various fluorescent dyes can be coupled to the 5'-end of the sense strand oligonucleotide to track transfection efficiency of the corresponding duplex.

* Larger synthesis scales are available on request.

siRNA

Custom siRNA Duplexes

Eurogentec has co-developed an exclusive siRNA design platform. PhD-level scientists of our design team use this reliable interface to design custom siRNA for any target of your choice.

Eurogentec **guarantees up to 80% minimum silencing** of your gene of interest with at least one of the 3 duplexes designed and synthesised. ■

Control siRNA Duplexes

In order to monitor your siRNA experiment conditions, Eurogentec provides siRNA control duplexes and kits including negative and positive controls necessary to validate your experiment.

■ **Negative controls** are siRNA molecules **presenting no homology with any known eukaryotic gene**. siRNA controls are already annealed and shipped in solution. The sequence is properly validated.

■ **Positive controls** consist of siRNA **directed against a range of endogenous and reporter genes**. They are available in 5nmol final quantities. Each control contains 1 siRNA duplex. All siRNA control duplexes are PAGE purified and 100% MALDI-TOF Mass Spectrometry controlled. The sequences are validated and published. ■

Catalogue control siRNA	
Name	Reference
Control siRNA duplex negative control	SR-CL000-005
Control siRNA duplex LaminB1 (human)	SR-CL001-005
Control siRNA duplex Vimentin (human)	SR-CL002-005
Control siRNA duplex NuMA (human)	SR-CL003-005
Control siRNA duplex Beta-actin (human)	SR-CL004-005
Control siRNA duplex Eg-5 (human)	SR-CL005-005
Control siRNA duplex Cdk-1 (human)	SR-CL006-005
Control siRNA duplex pGL2 luciferase (firefly)	SR-CL010-005
Control siRNA duplex pGL3 luciferase (firefly)	SR-CL011-005
Control siRNA duplex GFP (jellyfish)	SR-CL020-005

CUSTOM RNA OLIGO CAN BE ORDERED ONLINE VIA THE CUSTOM OLIGO CONFIGURATOR.

Principle of siRNA-mediated RNA interference. The annealed siRNA enter the cell (1). Once inside, double stranded RNA is recognised by the RISC complex. Sense strand siRNA is displaced and the mRNA anneal to the antisense siRNA fixed to the RISC complex (2). mRNA is digested (3) and the RISC complex containing the siRNA is then recycled to begin a new cycle (4).

miRNA

miRNA (for microRNA) are natural small non-coding RNAs forming short hairpins. They are implied in gene expression and RNA silencing.

Clear-MiR™ miRNA Inhibitors


Clear-MiR™ miRNA inhibitors are chemically modified antisense RNA oligonucleotides optimised to specifically **target miRNA molecule in cells**. ■

Add-MiR™ miRNA Mimics

Add-MiR™ oligonucleotides are custom double-stranded synthetic miRNA **mimicking the action of endogenous miRNAs**. ■

DID YOU KNOW

2'-O-Me RNA base and phosphorothioate links bring to the RNA oligo a higher stability and a resistance against nuclease. ■



>Note
Clear-MiR™ miRNA Inhibitors and Add-MiR™ miRNA Mimics are available with different labels and can be linked to cholesterol to increase cellular uptake. On request, peptides can also be covalently linked.

qPCR Probes

Eurogentec offers a wide range of fluorophores and quenchers in various combinations to fit any method and Real-Time thermocycler. ■

SPECIFICATIONS

Length: From 15 to 50 bases

Synthesis scale: 10 nmol • 40 nmol • 200 nmol • 1000 nmol • 2.5 µmol • 5 µmol • 10 µmol*

Backbone: DNA, LNA®, 2'-O-Me RNA and phosphodiester linkage

Modifications: 5': 6-FAM, HEX, Cy®3, TET, Cy5®...
3': TAMRA, DABCYL, BHQ™, DDQ...

Purifications: RP-HPLC or Dual HPLC

Quality Control: MALDI-TOF MS and analytical HPLC

Format: Dried

Packaging: 2 mL tube

Documentation: Technical data sheet

Probe Design: Available on request

Shipping: At room temperature

FOR COMMERCIAL USE

Access™ Dyes from Eurogentec are a class of high performance proprietary molecules with IP-friendly terms and conditions. ■

ALSO AVAILABLE

Eurogentec provides kits & reagents for **qPCR assays** (see Amplification brochure). ■

For **other modified oligonucleotides** please refer to the custom oligonucleotide chapter (p.9). ■

1 PROBE ORDERED = 2 PRIMERS OFFERED¹

>Note

1. Double-Dye or Molecular Beacon Probe. Each primer must be 15-30 DNA bases, unmodified and RP-Cartridge purified. The synthesis scale of the primers must be similar to the probe one (10, 40, 200 or 1000 nmol). Offer Valid only for orders placed on the Eurogentec e-commerce platform EOS. General conditions of sale will be applied.

* Larger synthesis scales are available on request




DOUBLE-DYE PROBES

EUROGENTEC OFFERS a large range of fluorescent dyes and quenchers including proprietary efficient molecules: HiLyte Fluor™ dyes and QXL™ quenchers.

LNA® Double-Dye probes

LNA® bases have a modification to the ribose backbone that locks the base in the C3'-endo position, which favors RNA A-type helix duplex geometry.

Compared to DNA Double-Dye probes, LNA®



Double-Dye probes exhibit higher thermal stabilities, specificity and reproducibility. They show better mismatch discrimination which allows the use of shorter probes.

Furthermore, LNA® offers the possibility to adjust Tm values of primers and probes in multiplex assays. ■




MOLECULAR BEACONS

EUROGENTEC IS a licensed supplier of Molecular Beacons and offers standard, wavelength-shifting and 2' O-Me RNA molecular beacon.

2' O-Methyl RNA Molecular Beacons

2' O-Methyl RNA probes perform better than DNA oligonucleotides. They are more nuclease resistant, have a higher affinity, specificity and hybridisation kinetics compared to DNA homologues.




Wavelength-Shifting Molecular Beacons

Wavelength-Shifting Molecular Beacons are brighter than standard Molecular Beacons because of an enhancement of the fluorophore signal. ■

ALSO AVAILABLE

Plexor™ primers. ■

LC HYBRIDISATION PROBES



Two LC hybridisation probes labelled with a single fluorescent molecule specifically recognise two adjacent sequences in the target DNA. When the probes are bound to the target sequence, the fluorescent signal is transferred from the donor to the acceptor, which starts to fluoresce. A 3' phosphate group is also added to prevent extension of the reporter probe by Taq DNA polymerase during the PCR cycles. ■

>Note

PNA FISH Probes are also available to detect chromosome aberrations in the centromer. Please see p.19 for more information.


MGB PROBES




EUROGENTEC PROVIDES high quality MGB probes perfectly suited for patient management [1]. MGB increases the Tm of a probe because of its minor groove binding ability. MGB probes are more specific, more efficient and more sensitive than standard double-dye probes.

We provide a complete offer with more than 15 dyes covering all qPCR channels - FAM,

- TET,
- AP5,
- Yakima Yellow®,
- Texas Red®,
- Cy®5,
- ROX,
- DragonFly™ Orange,
- ATTO
- HEX
- JOE



Our MGB Probes are RP-HPLC purified and can be delivered in 6, 20 and 50 nmol, dried or in solution (TE or H₂O). For maximal convenience, a 10 nmol dried aliquoted format is also available for the 20 and 50 nmol quantities at no additional cost. The probes are quality controlled by MALDI-TOF MS + HPLC and are available in IVD grade upon request. ■



The dispensing service, in line with ISO15189, brings to your in-house PCR/qPCR assays or commercial kits, a diagnostic-like grade fitting the highest quality standards. **Any size of routine assays to full kitting solutions** can be produced with a very high reliability, reproducibility and accuracy. This process saves set-up time and reduces reagent wastage, while keeping format flexibility.

High-throughput Dispensing service

Contact:
dispensing@eurogentec.com

Food
GMO
Clinical diagnostics
Veterinary pathogen
detection
Custom Buffer
Formulation

>Note

¹ Restriction of use in the following countries: US, CA, AU, CH, FR, UK, DE, IT, JP, SE, ES, CN. In these countries MGB probes must only be used for patient management. Use is free of limitation in other countries. End users are covered under Eurogentec's conveyed license for patient management.




FLUORESCENT DYES

QUENCHERS

Research & Commercial use

Research use only



ORDER
Oligonucleotides can be ordered via the web oligo configurator or on request at unique@eurogentec.com. For more information, please contact us at oligo@eurogentec.com

HIGH PERFORMANCE
Access™ Dyes and Quenchers* from Eurogentec are a class of high performance proprietary molecules with IP-friendly terms and conditions.

SEE ALSO
Compatibility cyler channel vs fluorescent dyes page 27.

USAGE
In bold blue: Dyes & Quenchers for commercial & research use
In black: Dyes & Quenchers for research use only (RUO)
In light grey: Not available at Eurogentec
Dark blue bands represent the indicative quenching range of recommended quenchers.
Gray bands represent the indicative quenching range of alternate quenchers.
Purple bars represent the optimal recommended quenching range.
Access™ Dyes & Quenchers

All wavelength values are in nm.

Note: Spectral properties of a dye or a quencher may vary when coupled to another molecule

BECAUSE YOUR EXPERIMENTS REQUIRE ALWAYS MORE CUSTOMISATION, UNIQUE OLIGONUCLEOTIDES BRING YOU THE PERFECT SOLUTION.

Unique™ oligos

⇒ HIGHLY COMPLEX OLIGOS

If you cannot find here the oligonucleotides that best suits your needs please think about unique oligos.

Eurogentec can synthesise highly complex oligonucleotides and can incorporate any of the modifications from recognised chemical suppliers (Glen research, TriLink BioTechnologies...).

Send us the specifications of your Unique™ oligonucleotides (sequence or length, chemistries, modifications, purifications, expected purity, synthesis scale or final amount, format, packaging...) or fill the form online (<https://secure.eurogentec.com/unique-oligonucleotides-quote-request-form.html>) and you will receive the corresponding information in terms of technical feasibility, pricing and turnaround times within 2 working days. ■

SPECIFICATIONS

Length: From 2 to 225 bases

Synthesis scale: Customised

Backbone: Usual or atypical chemistry

Modifications: Common or rare modifications

Purifications: SePOP desalting, RP-Cartridge·Gold™, HPLC, PAGE, Dual HPLC, UltraPureGold™

Quality Control: Adapted to your needs

Format: Adapted to your needs

Packaging: Adapted to your needs

Documentation: Technical data sheet custom documentation

Shipping: As defined by the customer

ALSO AVAILABLE

Custom Gene Synthesis

- From simple gene to highly complex sequence
- Up to 50 kbp
- Gene Optimisation
- 100% Guaranteed sequence
- Fast turnaround time

Contact: gene@eurogentec.com

More info in the [SmartGene Brochure](#)

Catalogue oligos

Name	Sequence	Bases	Tm (°C)	Ref.
16S rRNA For	AGA GTT TGA TCC TGG CTC AG	20	55.2	UN-PRO01-005
16S rRNA Rev	ACG GCT ACC TTG TTA CGA CTT	21	57.4	UN-PRO05-005
3' RACE PCR	GGC CAC GCG TCG ACT AGT AC	20	60.6	UN-PRO10-005
Anchored Oligo dT (20)	TTT TTT TTT TTT TTT TV	20	39.2	UN-PRO15-005
Anchored Oligo dT (22)	TTT TTT TTT TTT TTT TV N	22	42.8	UN-PRO20-005
Bluescript KS	TCG AGG TCG ACG GTA TC	17	53.3	UN-PRO25-005
Bluescript SK	CGC TCT AGA ACT AGT GGA TC	20	52.4	UN-PRO30-005
cDNA Cloning Primer	GGC CAC GCG TCG ACT AGT ACT TTT TTT TTT TTT TV	38	64.8	UN-PRO35-005
EGFP-C	CAT GGT CCT GCT GGA GTT CGT G	22	61.2	UN-PRO40-005
EGFP-N	CGT CGC CGT CCA GCT CGA CCA G	22	67.2	UN-PRO45-005
G3PDH For	ACC ACA GTC CAT GCC ATC AC	20	58.6	UN-PRO50-005
G3PDH Rev	TCC ACC ACC CTG TTG CTG TA	20	59.7	UN-PRO55-005
M13 Forward (-20)	GTA AAA CGA CGG CCA GT	17	53.0	UN-PRO60-005
M13 Forward (-41)	CGC CAG GGT TTT CCC AGT CAC GAC	24	65.5	UN-PRO65-005
M13 Reverse (-27)	CAG GAA ACA GCT ATG AC	17	47.3	UN-PRO70-005
M13 Reverse (-48)	AGC GGA TAA CAA TTT CAC ACA GG	23	57.2	UN-PRO75-005
Neomycin For	CTT GGG TGG AGA GGC TAT TC	20	55.6	UN-PRO80-005
Neomycin Rev	AGG TGA GAT GAC AGG AGA TC	20	54.0	UN-PRO85-005
Oligo dT, 15mer	TTT TTT TTT TTT TTT	15	29.7	UN-PRO90-005
Oligo dT, 16mer	TTT TTT TTT TTT TTT T	16	32.1	UN-PRO95-005
Oligo dT, 18mer	TTT TTT TTT TTT TTT TTT	18	36.0	UN-PR100-005
Oligo dT, 20mer	TTT TTT TTT TTT TTT TTT TT	20	39.1	UN-PR105-005
PCMV Forward	CGC AAA TGG GCG GTA GGC GTG	21	64.8	UN-PR110-005
pET 3'	CTA GTT ATT GCT CAG CGG	18	50.6	UN-PR115-005
pET 5' (T7)	TAA TAC GAC TCA CTA TAG G	19	45.3	UN-PR120-005
pET Upstream	ATG CGT CCG GCG TAG A	16	56.7	UN-PR125-005
pGEX 3'	CCG GGA GCT GCA TGT GTC AGA GG	23	65.2	UN-PR130-005
pGEX 5'	GGG CTG GCA AGC CAC GTT TGG TG	23	67.0	UN-PR135-005
ROSA26 Promoter For	AAA GTC GCT CTG AGT TGT TAT	21	53.2	UN-PR140-005
ROSA26 Promoter Rev	GGA GCG GGA GAA ATG GAT ATG	21	56.3	UN-PR145-005
SP6 Promoter	TAC GAT TTA GGT GAC ACT ATA G	22	50.0	UN-PR150-005
SP6 Upstream	ATT TAG GTG ACA CTA TAG	18	42.8	UN-PR155-005
T3 Promoter	AAT TAA CCC TCA CTA AAG GG	20	50.4	UN-PR160-005
T7 Promoter	TAA TAC GAC TCA CTA TAG GG	20	48.3	UN-PR165-005
T7 Terminator	GCT AGT TAT TGC TCA GCG G	19	54.1	UN-PR170-005

⇒ UNIVERSAL PRIMERS

Universal primers are complementary to nucleotide sequences that occur very commonly in specific sets of DNA molecules and cloning vectors. Thus, they are able to **bind to a wide variety of DNA templates.** ■

SPECIFICATIONS

Quantity: 1 OD/5 nmol

Backbone: DNA

Modifications: None

Purifications: RP-HPLC

Quality Control: MALDI-TOF MS + CGE

Format: Dried

Packaging: 2 mL tube

Documentation: Technical data sheet

Shipping: At room temperature

CONTACT US

unique@eurogentec.com

➔ PNA FISH

In principle, fluorescence *in situ* hybridisation (FISH) should be able to provide information on the telomere length of individual chromosomes. The efficiency of conventional labelled oligos is not sufficient to be extended beyond qualitative studies of TTAGGG repetitions in chromosomes of various species. PNA chemical structure brings a **higher sequence specificity**, an improved **stability**, better **reproducibility**, and lower background noise. Due to the higher T_m of PNA/DNA duplexes, short (18-mer) telomere PNA (CCCTAA)₃ are now widely used. ■

SPECIFICATIONS

Length: 18 bases
Quantity: 5 nmol
Backbone: PNA
Modifications: FAM • Cy3[®] • Cy5[®] • FITC • TMR

➔ CALIBRATION OLIGOS

Dye-labelled calibration oligos are a set of 5' fluorescent dT10 oligonucleotides recommended to calibrate some real-time qPCR thermocyclers. Calibration is a preliminary step indicated to adjust fluorescent signal analysis. qPCR Dye Calibration Oligos' enables the thermocycler to recognise the spectra of each single dye and to control signal overlap that may occur in multiplexed assays particularly. ■

SPECIFICATIONS

Length: 10 bases
Quantity: 5 nmol
Backbone: DNA
Modifications: AP5, Yakima Yellow[®], HEX, Dragonfly Orange[™], TET, JOE, HiLyte[™] Fluor 647, ROX
Purification: RP-HPLC
Quality control: MALDI-TOF MS
Format: Dried
Packaging: 2 mL tube
Documentation: Technical data sheet
Shipping: At room temperature

ORDERING INFORMATION					
Name	Quantity	#Cat	Name	Quantity	#Cat
C-Rich Telomere Probes			Centromere Probes		
TelC-FAM	5 nmole	PN-TC001-005	Cent-Cy3	5 nmole	PN-CN050-005
TelC-Cy3	5 nmole	PN-TC050-005	Cent-FAM	5 nmole	PN-CN001-005
TelC-Cy5	5 nmole	PN-TC055-005	Cent-Cy5	5 nmole	PN-CN055-005
TelC-Alexa488	5 nmole	PN-TC060-005	Cent-Alexa488	5 nmole	PN-CN060-005
TelC-FITC	5 nmole	PN-TC011-005	Cent-FITC	5 nmole	PN-CN011-005
TelC-TAMRA	5 nmole	PN-TC030-005	Cent-TAMRA	5 nmole	PN-CN030-005
TelC-Alexa647	5 nmole	PN-TC020-005	Cent-Alexa647	5 nmole	PN-CN020-005
TelC-Biotin	5 nmole	PN-TC040-005	Cent-Biotin	5 nmole	PN-CN040-005
G-Rich Telomere Probes			Centromere Protein B Probes		
TelG-FAM	5 nmole	PN-TG001-005	CENPB-FAM	5 nmole	PN-CP030-005
TelG-Cy3	5 nmole	PN-TG050-005	CENPB-Cy3	5 nmole	PN-CP050-005
TelG-Cy5	5 nmole	PN-TG055-005	CENPB-Cy5	5 nmole	PN-CP055-005
TelG-Alexa488	5 nmole	PN-TG060-005	CENPB-Alexa488	5 nmole	PN-CP060-005
TelG-FITC	5 nmole	PN-TG011-005	CENPB-FITC	5 nmole	PN-CP011-005
TelG-TAMRA	5 nmole	PN-TG030-005	CENPB-TAMRA	5 nmole	PN-CP001-005
TelG-Alexa647	5 nmole	PN-TG020-005	CENPB-Alexa647	5 nmole	PN-CP020-005
TelG-Biotin	5 nmole	PN-TG040-005	CENPB-Biotin	5 nmole	PN-CP040-005

ORDERING INFORMATION					
Name	Sequence	Modification 5'	Bases	Abs/Em (nm)	Reference
AP5-T10	TTTTTTTTTT	AP5	10	527/549	UN-CT001-005
YY-T10	TTTTTTTTTT	YY	10	530/549	UN-CT005-005
HEX-T10	TTTTTTTTTT	HEX	10	535/556	UN-CT010-005
DFO-T10	TTTTTTTTTT	DFO	10	554/576	UN-CT015-005
TET-T10	TTTTTTTTTT	TET	10	521/536	UN-CT020-005
JOE-T10	TTTTTTTTTT	JOE	10	529/555	UN-CT025-005
HL647-T10	TTTTTTTTTT	HL647	10	650/675	UN-CT030-005
ROX-T10	TTTTTTTTTT	ROX	10	575/602	UN-CT035-005

Glen Synthesis reagents

Eurogentec is an **authorised distributor for Glen research** products in most of European countries. We distribute a large range of reagents for the DNA/RNA oligonucleotide synthesis including phosphoramidites, solid supports for oligonucleotide purifications, labelling dyes and modifications. ■

- CE Phosphoramidites, synthesis columns and solvents/reagents
- Other monomers
- Minor bases
- Modification and labelling
- RNA Synthesis
- And many more

Additional Services

Additional QC

MALDI-TOF Mass Spectrometry: This method provides the most precise information about the length, deprotection-product and the presence of labels for modified oligonucleotides over a broad range of lengths (up to 60 bases).

RP-UHPLC: This is a very efficient technique giving quantitative information about the purity level of oligonucleotides from 15 to 40 bases long.

IEX-UHPLC: This technique is particularly adapted to quantify the purity level of oligonucleotides from 15 to 40 bases long.

Capillary Gel Electrophoresis (CGE): This method is adapted to assess very precisely the purity of oligonucleotides longer than 40 bases (on request).

Fluorescence analysis: This non-destructive physical technique provides qualitative information about your fluorescent oligonucleotides. ■

Format

Dried: All the synthesised oligonucleotides are dried by default (except SePOP unmodified oligonucleotides from 15 to 39 bases).

In solution: You may select the nature of the reconstitution buffer (H₂O or TE), the volume of the reconstitution buffer (from 50 to 1000 µl) or/and the final oligonucleotides concentration (from 5 to 250 µM).

Annealed: siRNA or cloning linkers are annealed by default.

Mixed: Similar amounts of forward and reverse oligonucleotides can be mixed in a single tube. ■

Packaging

2 mL tube: By default, each oligonucleotide is provided in individual 2 mL tube. Higher volume can be delivered on request (15 mL, 50 mL)

96-well plates: Cluster tubes, well plate and deep well plate are available.

384-well plates: Specially suitable for high throughput experiments requiring more than 96 oligonucleotides.

Aliquoting: All the oligonucleotides in solution can be split in small aliquots of desired volume (from 50 to 1000 µl). ■

Shipping

Your oligonucleotides can be express shipped in 24 hours upon request (see page 23 for more details). ■

Design

We continuously update our software and design rules to reflect the latest scientific developments as well as integrate customer requirements. This service includes primers, Double-Dye Oligonucleotides, Molecular Beacons, siRNA design, miRNA inhibitors... ■

Increase the quality of your oligo for demanding applications

Hospital and commercial kits

Track™ oligonucleotides

Track™ oligonucleotides offer a higher traceability (and other quality assets) in the production process than life science research oligonucleotides.

cGMP oligonucleotides

cGMP oligonucleotides ensure exceptional product quality by manufacturing in classified cleanrooms and use of an ISO 13485-certified and GMP compliant QMS.

Therapeutics fields

Pre-clinical oligonucleotides

Large scale pre-clinical oligonucleotides are manufactured in cleanrooms and delivered with appropriate documentation. Additional QC tests such as endotoxin level are offered. For more information download our IVD brochure.

www.eurogentec.com/invitrodiagnostics.html

OLIGO GRADE						
Lab Services Hospital & Clinical Labs		Discovery	Routine Assays	Contamination Sensitive Assays		
Diagnostic Companies		Discovery	Feasibility	Prototyping	Validation	Commercialisation
Oligonucleotide Grades	Research	Track	Pre-Diagnostic	Diagnostic		
Process						
Dedicated Account Contact Person	Option	✓	✓	✓	✓	
Customised Fill & Finish	Option	Option	✓	✓	✓	
Quality Management						
ISO 9001 Certification	✓	✓	✓	✓	✓	
ISO 13485 Certification	-	-	✓	✓	✓	
Qualification/Validation [Equipment & Method]	-	Partial	✓	✓	✓	
Control						
Quantification	Single	Dual	Triple	Triple		
Stringent QC Tests (validated)	-	✓	✓	✓		
Traceability	Partial	Documented	Documented	Full documented		
Batch Record [Archived for 5 years]	-	-	Partial	Full		
Classified Cleanroom	-	-	✓	✓		
Certificate of Analysis [CoA]	-	✓	✓	✓		

ANNEXES

SYNTHESIS SCALE VS GUARANTEED YIELD

The **synthesis scale** refers to the **amount of raw material** used to start the synthesis of oligonucleotides. The **yield** corresponds to the amount of **final product**

recovered at the end of the synthesis and purification processes. The length, the sequence, the type/number of modifications and the purification, strongly

influence the reaction yield. Based on that, Eurogentec defined a minimum guaranteed yield in nmoles for all product categories (see table below). The minimum guaranteed yields

represent only a reference because the delivered quantities may vary.

GOOD TO KNOW

All allowed purifications are represented in this table. To select the recommended purification according to your applications and modifications, please refer to p.8.

		Synthesis scale (nmol)																																
		10	40	200	1000	2500	5000	10000	20000																									
		Purification																																
Range	Product	Length	SePOP	RP-Cartridge-Gold™ HPLC (RP or EX)	SePOP	RP-Cartridge-Gold™ HPLC (RP or EX)	PAGE ¹³	Dual HPLC	UltraPureGold™	SePOP	RP-Cartridge-Gold™ HPLC (RP or EX)	PAGE ¹³	Dual HPLC	UltraPureGold™	SePOP	HPLC (RP or EX)	Dual HPLC	SePOP	HPLC (RP or EX)	Dual HPLC	SePOP	HPLC (RP or EX)	Dual HPLC											
Custom Oligonucleotides	Non-Modified (DNA only)	5-9	-	-	-	-	-	-	-	60	50	30	20	15	-	180	100	80	40	40	-	450	200	100	900	400	200	1800	800	400	-	-	-	
		10-19	-	-	-	-	-	-	3	-	70	60	45	30	23	15	200	140	100	70	50	30	500	250	125	1000	500	250	2000	1000	500	4200	2100	1050
		20-39	5	4	-	20	16	10	4	2	60	50	30	20	15	-	190	120	90	40	45	-	475	225	115	1000	500	250	2000	1000	500	4200	2100	1050
		40-59	3	2	-	10	8	5	2	1	30	25	15	12	7	-	115	60	45	20	20	-	285	110	55	600	230	115	1200	460	230	2500	1000	500
		60-79	-	-	-	-	6	-	2	-	18	-	8	-	-	-	40	-	14	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
		80-99	-	-	-	-	-	-	1	-	-	-	-	3	-	-	-	-	5	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
100-139	-	-	-	-	-	-	-	-	-	-	-	2	-	-	-	-	3	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
Real-Time qPCR Probes	Modified ⁽¹⁾ (including DNA, RNA, 2' O-Me RNA, LNA and phosphorothioate linkages)	5-9	-	-	-	-	-	-	-	-	12	-	-	-	-	25	-	12	-	-	-	60	30	-	125	60	-	250	125	-	-	-		
		10-19	-	-	-	12	6	5	4	1	35	20	17	15	8	-	70	40	35	30	15	-	90	45	-	190	95	-	380	190	-	-	-	
		20-59	-	-	-	8	5	4	3	1	20	15	12	10	6	-	45	35	25	20	12	-	65	30	-	135	65	-	275	130	-	600	275	
		60-139	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
Real-Time qPCR Probes	Double-Dye probes ⁽²⁾	8-38	-	-	<2 ⁽⁴⁾	-	4	-	-	-	12	-	-	-	-	25	-	-	-	-	-	65	-	-	135	-	-	275	-	-	600	-		
		Molecular Beacons	32-50	-	-	1	-	-	-	-	4	-	-	-	-	-	12	-	-	-	-	-	30	-	-	65	-	-	130	-	-	275	-	
		MGB Taqman Probes	8-30	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
Delivered quantity: 6, 20 or 50 nmol																																		
RNAi Oligonucleotides	siRNA Duplexes Non-Modified ⁽⁵⁾	21-27	7	-	3	22	-	12	-	60	-	40	-	-	200	-	80	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
	siRNA Duplexes Modified ⁽¹⁾	21-27	7	-	3	22	-	12	-	60	-	40	-	-	200	-	80	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
NGS Oligonucleotides	RP-Cartridge purified	20-85	Minimum delivered quantity: 10 nmol																															
	RP-HPLC purified																																	
Universal Primers	-	15-38	-	-	-	-	5	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
Unique Oligonucleotides	-	2-225	On request - please contact us at unique@eurogentec.com																															

Post-synthesis modifications may yield 50% less than the above stated values.

Table: (1) Between 5 and 59 bases length single-modified Oligonucleotides. Eurogentec does not provide minimum guaranteed yield for modified oligonucleotides greater than 59 bases. Post-synthesis modifications are not compatible with SePOP and RP-Cartridge-Gold™ purification. A lower yield may result from poly-modifications and/or strong secondary structures.

(2) Double-Dye probes only result from the combination of a 5' fluorescent dye and a 3' quencher.

(3) Except for oligonucleotides with GC-rich regions.

(4) Only available for Double-Dye FAM-TAMRA 10 nmol and FAM-BHQ1 10 nmol.


(5) Non-modified siRNA's only include 3' dTdT overhang.

List of the post-synthesis modifications

- > 5' Alexa Fluor* (350, 430, 488, 500, 514, 532, 546, 555, 568, 594, 610, 633, 647, 660, 680, 700 and 750)
- > 5' ATTO (390, 425, 465, 488, 495, 520, 532, 550, 565, 590, 594, 610, 620, 633, 635, 647N, 655, 680, 700, 725 and 740)
- > 5' BODIPY* (530/550, FL and TR)
- > 3', 5' and dT Cascade Blue*
- > 3' and dT Cy* (3, 3.5, 5 and 5.5)
- > 3', 5', dR and dT Digoxigenin
- > 5' Dragonfly Orange*
- > 5' DY- (681, 781 and 782)
- > dR 6-FAM

- > dR and dT HEX
- > 5' HiLyte™ Fluor (405, 488, 555, 594, 647, 680 and 750)
- > 3', dR and dT JOE
- > 5' Marina Blue*
- > 5' Oregon Green* (488 and 488 X)
- > 5' Pacific Blue*
- > 3' QXL*
- > 3', 5', dR and dT Rhodamine 6G
- > 3', 5', dR and dT ROX
- > 5' TAMRA
- > dR and dT TET
- > 3', 5', dR and dT Texas Red*

For more information, please contact our Oligo Centre at: oligo@eurogentec.com or visit our website: www.eurogentec.com



FOR PCR DETERMINE THE RIGHT SYNTHESIS SCALE

Final Oligo Concentration	50 nM	150 nM	300 nM	600 nM	900 nM	Minimum quantity to order*
Average number of Reactions (total volume 100µL)						
100	30	15	8	6		0.5 nmol
1000	300	165	80	55		5 nmol
5000	1650	830	415	275		25 nmol
10 000	3300	1660	830	555		50 nmol
100 000	33 300	16 660	8330	5555		500 nmol

*Please select in the minimum guaranteed yield table the synthesis scale corresponding to the desired minimum quantity.

CHECK YOUR SHIPPING METHODS

The delivery time depends on the specifications of your oligonucleotides (see table below).

Eco-Logik Delivery

■ By local Mail to reduce the global ecological impact. Receive your oligonucleotides in your mailbox. Available for Belgium, France and Monaco

Express Delivery

■ All oligonucleotides
■ By Express courier to receive your oligonucleotides as fast as possible (24 to 48 hours) in your hands.

■ Same day shipping option
- For orders received before 10.00 AM (Central European Time)
- For custom oligonucleotides (max 24), 10/40 nmol scale,
5-30 DNA bases, unmodified, SePOP desalted or RP-Cartridge purified.

Delivery times (in working day)

Range	Product	Length	Purification					
			SePOP	RP-CartridgeGold™	HPLC (RP or IEX)	PAGE	Dual HPLC	UltraPureGold™
Custom Oligonucleotides	Non-Modified (DNA Only)	5-9	-	4-5	5	6	7	7
		10-39	2-3	4-5	5	6	7	7
		40-59	5	6	7	8	-	9
		60-79	-	6	-	8	-	9
		80-139	-	-	-	10	-	11
	Modified (including DNA, RNA, 2' O-Me RNA & LNA)	10-39	5	7	7	8	9	9
Real-Time qPCR Probes	Double-Dye Probes	8-38	-	-	7	-	-	-
	Molecular Beacons	32-50	-	-	12-15	-	-	-
	MGB Taqman Probes	8-30	-	-	5-7	-	-	-
RNAi Oligonucleotides	siRNA Duplexes	21-27	5-7	-	9-10	10-11	-	-
NGS Oligonucleotides	-	20-85	-	4-6	5-7	-	-	-
Universal Primers	-	15-38	-	-	2-3	-	-	-
Unique Oligonucleotides	-	2-225	On Request					

For large order or Unique Oligonucleotides, please feel free to contact us at oligo@eurogentec.com to receive more details in terms of delivery schedules. 5'AP, BSA, HRP or SBP Conjugation : 3-5 WD Extra. Additional Purification or Services: 2 WD Extra ; Fax Ordering: 1 WD Extra

RECEIVE YOUR DOCUMENTATION

Each oligonucleotide is provided with a technical data sheet. Other documentations could be added depending on the oligonucleotide type. All the documents are sent as pdf files to your shipping email address.

	TDS	MS ⁽¹⁾	UHPLC	CGE ⁽³⁾
Custom Oligonucleotides	Unmodified	✓		
	Modified	✓	✓ ⁽⁴⁾	
	UltraPureGold™	✓	✓	✓
Real-Time qPCR Probes	✓	✓	✓	
RNAi Oligonucleotides	✓	✓ ⁽⁴⁾		
Universal Primers	✓	✓		
Unique Oligonucleotides	✓	✓	✓ ⁽⁵⁾	✓ ⁽⁵⁾
NGS oligos	✓	✓		
Calibration oligos	✓	✓		

TDS: Technical Data Sheet; MS: Mass Spectrometry; HPLC: High Performance Liquid Chromatography; Ultra Performance Liquid Chromatography; CGE: Capillary Gel Electrophoresis.
(1) Always provided up to 60 bases long Oligonucleotides.
(2) If applicable.
(3) Can be substituted by another analytical QC
(4) Except for SePOP desalted oligonucleotides.
(5) Optional.
For technical reasons this general rule may be adapted to provide you with the most suitable and useful documentation.

HOW TO STORE YOUR OLIGO

Handling information

Products	Format	Storage	Stability**
Custom Oligonucleotides	Dried TE Buffer (pH 8) or dH ₂ O	RT -20 °C	18 months 24 months
Real-Time qPCR Probes	Dried TE Buffer (pH 8)* or dH ₂ O	RT -20 °C	18 months 24 months
RNAi Oligonucleotides	Dried RNase-free Buffer (pH 7.5)	RT -20 °C	18 months 24 months
Catalogue Primers PNA FISH Probes / Custom PNA	Dried	RT	18 months

* Except for Cy⁵ dye labelled oligonucleotides (pH7)

** Please protect from light and avoid freeze/thaw cycles.

Please note that depending on sequence and modifications, the stability of the oligos may vary substantially versus the values given above, which should therefore be considered as indicative.

RECONSTITUTE YOUR OLIGO

1. Spin the tube briefly to collect the pellet in the bottom of the tube.
2. Add an appropriate volume of recommended buffer.
3. Allow the tube to stand a few minutes.
4. Vortex the tube for 15 seconds and spin briefly.
5. Refer to the dedicated technical data sheet for more information.



QUANTIFY YOUR OLIGO

To quantify your oligonucleotides, make an aliquot of the resuspended oligonucleotides to a final volume of 1 mL of dH₂O and vortex for a few seconds. Measure the absorbance of this dilution at 260 nm (A₂₆₀). Use the formula below to calculate the concentration of oligonucleotides in your stock solution. This formula is valid for an absorption of A₂₆₀ ≤ 1.2.

Concentration in µg/mL = A₂₆₀ × dilution factor × Weight per OD of stock solution (in µg / OD).

1 OD₂₆₀ (Optical Density) unit is defined as the amount of oligonucleotide which, when dissolved in a volume of 1.0 mL, results in an absorbance of 1.0 when measured at 260 nm in a 1 cm path-length quartz cuvette. 1 OD₂₆₀ unit corresponds to approximately 33 µg of single strand DNA. These relationships, however, can be inaccurate for short fragments of DNA, such as oligonucleotides. Base composition and even linear sequence will affect optical absorbance. Hence the precise value of the OD to mass relationship is unique for each oligonucleotide.

MEASURE

1.0 OD₂₆₀ of CCCCCCCCCC (10 bases) equals 39 µg
whereas 1.0 OD₂₆₀ of AAAAAAAAAA (10 bases) equals only 20 µg.

We carefully measure the OD value for your custom oligonucleotide by measuring the absorption at 260 nm using UV spectrophotometer. This information is provided on the oligonucleotide technical data sheet as the number of OD₂₆₀ units. The amount of oligonucleotide expressed in nanomoles and micrograms is derived from the OD measurement.

CALCULATE

Calculate the number of nanomoles present given an OD reading and extinction coefficient:

$$\text{Nanomoles} = (\text{OD}_{260} / \epsilon_{260}) \times 10^6$$

Example:

1 OD₂₆₀ unit of primer M13 Forward,
5'-GTA AAA CGA CGG CCA GTG-3'
Molar extinction coefficient (ε₂₆₀) = 182.800 L / (mole x cm)
Nanomoles = (1.0 / 182.800) × 10⁶ = 5.47 nmoles

CONVERT

Convert the amount in nanomoles to micrograms:

$$\text{Micrograms} = \text{Molecular Weight} \times \text{Nanomoles} \times 10^{-3}$$

Example:

1 OD₂₆₀ unit of primer M13 Forward,
5'-GTA AAA CGA CGG CCA GTG-3'
Molecular Weight = 5558.7
Micrograms = 5558.7 × 5.47 × 10⁻³ = 30.4 µg

CALCULATE THE DATA

THE MOLAR EXTINCTION COEF.

$$\epsilon_{260} = 2 \times \left(\sum_1^{n-1} \epsilon_{\text{Nearest Neighbour}} \right) - \sum_2^{n-1} \epsilon_{\text{Individual}} + \sum_1^n \epsilon_{\text{Modification}}$$

where $\epsilon_{\text{Nearest Neighbour}}$ is the nearest neighbour constant for a pair of bases, $\epsilon_{\text{Individual}}$ is the constant for an individual base, and n is the length of the oligonucleotide.

THE MOLECULAR WEIGHT

$$\text{Anhydrous MW (g/mol)} = \sum_{\text{Individual Base}} \text{MW} + \sum_{\text{Individual Mods}} \text{MW} - 63.98 + 2.016$$

For DNA bases:

MW dA = 313.21; MW dC = 289.18; MW dG = 329.21; MW dT = 304.20; MW dU = 290.17; MW dl = 314.19

For RNA bases:

MW DNA counterpart + 16.
When determining the weight of Uracil (rU) start with dU and not dT


For LNA bases:

MW DNA counterpart + 16 (+42 for dC)

For 2' O-Methyl bases:

MW DNA counterpart + 30.03.
When determining the weight of mU start with dU and not dT

For phosphorothioated bases: MW DNA counterpart + 16.06



WRITE YOUR SEQUENCE IN YOUR ORDER

IUB CODE

ACGTA = DNA
(ACGUA) = RNA
[ACGUA] = 2' O-Me RNA
<ACGTA> = 2' O-MOE
{ACGTA} = LNA®
A*C*(G*U*A) = Phosphorothioate links

Mixed bases (also known as degenerate or wobble bases) follow the IUB codes:

D=A/G/T
M=A/C
H=A/C/T
I = Inosine = Universal base
W=A/T
R=A/G
Y=C/T
V=A/C/G
S=C/G
K=G/T
N=A/G/C/T
B=C/G/T

Oligonucleotide synthesised with mixed bases gives a final product that is a heterogeneous population of distinct species. MW, Tm and extinction coefficient may be strongly affected by mixed base addition. Rather than reporting the various values for each component, a single value is given.



>Note

ε modification is not known for all modifications.

COMPATIBILITY CYCLER CHANNELS VS FLUORESCENT DYES

Thermocycler	Blue channel	Channel 1	Channel 2	Channel 3	Channel 4	Channel 5	Channel 6	Channel 7
GeneAmp SDS 5700		FAM	AP5/YY/DE/VC/TET*	DF0/TAMRA/NEO	ROX			
ABI Prism SDS 7000		FAM	AP5/YY/DE/VC/TET*	DF0/TAMRA/NEO	ROX			
ABI Prism SDS 7700		FAM	AP5/YY/DE/VC/TET*	DF0/TAMRA/NEO	ROX			
ABI Prism SDS 7900 HT		FAM	AP5/YY/DE/VC/TET*	DF0/TAMRA/NEO	ROX			
ABI Prism SDS 7300		FAM	AP5/YY/DE/VC/TET*	DF0/TAMRA/NEO	ROX			
ABI Prism SDS 7500		FAM	AP5/YY/HEX/HEX/DE/VC/TET*	DF0/TAMRA/Cy3/NEO	ROX/TR	Cy5/HL647		
QuantaStudio 3.5.6.7		FAM	AP5/YY/HEX/HEX/DE/VC/TET*	DF0/TAMRA/Cy3/NEO	ROX/TR	Cy5/HL647		
QuantaStudio 7		FAM	AP5/YY/HEX/HEX/DE/VC/TET*	DF0/TAMRA/Cy3/NEO	ROX/TR	Cy5/HL647		
VIAT		FAM	AP5/YY/HEX/HEX/DE/VC/TET*	DF0/TAMRA/NEO	ROX/TR*	LIZ/ATTO 633		
StepOne		FAM	AP5/YY/HEX/HEX/DE/VC/TET*	ROX				
StepOnePlus		FAM	AP5/YY/HEX/HEX/DE/VC/TET*	DF0/TAMRA/Cy3/NEO	ROX			
CyclerQ		FAM	AP5/YY/HEX/Cy3/TET	DF0/TAMRA/Cy3/NEO	ROX/TR	Cy5		
MyQ		FAM		DF0/TAMRA/Cy3/NEO	ROX/TR	Cy5		
Q5		FAM	AP5/YY/HEX/DE/VC/TET	DF0/TAMRA/Cy3/NEO	ROX/TR	Cy5		
CFX 96		FAM	AP5/YY/HEX/DE/VC/TET	ROX/TR	Cy5	Cy5.5		
MiniOpticon		FAM	AP5/YY/HEX/TET*					
DNA Engine Opticon1		FAM						
DNA Engine Opticon2		FAM	AP5/YY/HEX/TAMRA/VC/TET					
Chrome 4		FAM	AP5/YY/HEX/TAMRA/Cy3/DE/VC/TET	ROX/TR	Cy5			
MX300P (choice of 4 filters)		FAM	TET	YY/HEX/DE/VC	Cy3	DF0/TAMRA/NEO	TR/ROX	Cy5
MX300EP (choice of 5 filters)		FAM	TET	YY/HEX/DE/VC	Cy3	DF0/TAMRA/NEO	TR/ROX	Cy5
MX400T (choice of 4 filters)		FAM	TET	YY/HEX/DE/VC	Cy3	DF0/TAMRA/NEO	TR/ROX	Cy5
AriaMk (choice of up to 6 filters)		FAM	HEX	Cy3	ROX	Cy5		
Mastercycler ep realplex2		FAM	AP5/YY/HEX/DE/VC/TET					
Mastercycler ep realplex4		FAM	AP5/YY/HEX/DE/VC/TET	DF0/TAMRA	ROX			
LightCycler 1.5		FAM	AP5/HEX/YY/DE/VC/TET	Cy5				
LightCycler 2		FAM	AP5/HEX/YY/DE/VC/TET	TR/LC Red 600	ATTO 620/LC Red 640	Cy5/LC Red 670	LC Red 705/ATTO 680	
LightCycler 480 I		FAM	AP5/YY/HEX/DE/VC	ROX/TR/LC Red 610	ATTO 620/LC Red 640	Cy5		
LightCycler 480 II		FAM	AP5/YY/HEX/DE/VC	ROX/TR/LC Red 610	ATTO 620/LC Red 640	Cy5/Cy5.5		
LightCycler 96		FAM	AP5/YY/HEX/VC	TR/LC Red 610	Cy5			
SmartCycler1		FAM	AP5/YY/Cy3/DE/VC/TET	TR	Cy5			
SmartCycler2		FAM	AP5/Cy3/YY/DE/VC/TET	TR/ROX	Cy5			
Rotor-Gene 2000 / 3000		FAM	AP5/YY/DE/VC/TET	TAMRA/ROX/Cy3.5/TR	Cy5			
Rotor-Gene 6000		FAM	AP5/YY/HEX/DE/VC/TET	TAMRA/ROX/Cy3.5/TR	Cy5	ATTO 680		
PKiReal		FAM	AP5/HEX/YY	ROX/TR	Cy5			
Qower		FAM	AP5/YY/HEX/DE/VC	DF0/TAMRA/NEO	ROX/TR/Cy3.5	Cy5	Cy5.5	

In grey = Not available at Eurogentec

In blue = Recommended by Eurogentec

YY = Yakima Yellow
DF0 = Dragonfly Orange

* perform a dye calibration for optimal results
For complementary information, please refer to instrument manufacturer technical guide or contact us at scientific-support@eurogentec.com

HOW TO ORDER ON LINE

WWW.EUROGENTEC.COM

SINGLE ORDER

1. Connect to www.eurogentec.com
2. Click on the oligonucleotide tab of the order centre screen
3. Select the oligonucleotide type (Custom, Probes, RNAi...)
4. Fill the configurator with your oligonucleotide specifications
5. Add your oligonucleotide into your cart and finalise your order

MULTIPLE/ BATCH ORDERS

1. Connect to www.eurogentec.com
2. Click on the oligonucleotide tab of the order centre screen
3. Select the Multiple/Batch Order
4. Download the Excel File and fill it in
5. Upload the completed file on the Eurogentec website: www.eurogentec.com

HOW TO PAY

POSTPAID SYSTEM

One order / one invoice
You place an order of 1 or multiple oligonucleotides and you receive the invoice corresponding to this order.

PREPAID SYSTEM Oli&GO™

One invoice for multiple orders
You place a defined amount on your Oli&GO™ account. You receive an invoice corresponding to this amount. You can use this amount over time.



ECONOMIC Exclusive Oli&GO™ prices.

EASY Only one invoice for multiple oligo's orders spread over time.

CONVENIENT Oligo orders scheduled and tracked on line.

SECURED One administrator can give restricted access to multiple users.

With the realtime integrated counter, you keep an eye on your budget.

Scheduling your oligonucleotide orders allows reducing the number of parcels sent and decreases your shipping costs.

	Administrator Full of privileges to control the system	User Restricted access
Functionalities		
Use 1 or more Accounts	✓	✓
Buy Oligonucleotides	✓	✓
Receive an Order Confirmation	✓	✓
Receive the Related Documentation	✓	✓
Rename the Account	✓	✗
Add/Remove User(s)*	✓	✗
Define/Update Shipping Address	✓	✗
Reload the Account(s)	✓	✗
Schedule Orders (Day/Time)	✓	✗

*Multiple users can be defined per account

SHIPMENT GROUP

HOW TO REDUCE MY SHIPPING FEES

With the shipment group option, all the labs from the same institution can group their shipments to benefit from free (or reduced) shipping cost.

ECONOMIC
Free shipping of your oligos¹.

ECOLOGIC
Reduction of the number of parcels sent.

FAST
Shipping of the oligos as soon as they are ready².

RELATED PRODUCTS

Custom genes		
dNTPs		
dNTP Mix	1x 20 µmoles	NU-0010-10
dNTP Set	4 x 25 µmoles each	NU-0020-50
Takyon™ qPCR kits		
Test your free sample, visit www.eurogentec.com/qpcr-takyon.html		
DNA purification kits (100 preps)		
SMARTPure PCR Kit		SK-PCPU-100
SMARTPure Gel Kit		SK-GEPU-100
SMARTPure Plasmid Kit		SK-PLPU-100
DNA extraction kit (100 preps)		
SMARTExtract DNA kit		SK-DNEX-100
Agarose		
Agarose - 500g		EP-0010-05
Agarose small fragment (125g)		EP-0020-10
AgaTabs - 300 tablets		EP-0030-15
MW markers		
SmartLadder	(200 to 10000 bp)	MW-1700-10
SmartLadder SF	(100 to 1000 bp)	MW-1800-04
Electrophoresis devices		
Mupid®-One	EU cable	MU-0041-
	UK cable	MU-0041+
SmartViewer for Mupid®		MU-0101
SmartIlluminator		MU-0201

>Note The online tracking allows you to check the statements of your oligonucleotide orders at any time.

>Note ¹Minimum order amount required. ²Depending on the order quantity, we can determine a delivery plan. See our detailed shipping conditions on <http://www.eurogentec.com/shipping-conditions.html>

- For Roche Enzyme Products Conveying Polymerase Patent Rights, Chemically-Modified Hot-Start Polymerase Reagents and Kits, Licensed dsDNA-Binding Dye Research Kits and Licensed 5' Nuclease Research Kits
 - Use of this product is covered by applicable US patents and corresponding patent claims outside the US held by Roche. The purchase of this product includes a limited, non-transferable immunity from suit under such patent claims for using only this amount of product for the purchaser's own internal research. No right under any other patent claim, no right to perform any patented method, and no right to perform commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is conveyed expressly, by implication, or by estoppel. This product is for research use only. Additional rights for use in non-diagnostic fields require a license from Life Technologies. Diagnostic uses require a separate license from Roche. Further information on purchasing licenses may be obtained by contacting Out-licensing at Thermo Fisher Scientific, 5791 Van Allen Way, Carlsbad California 92008, outlicensing@lifetech.com.
- For products subject to Section 2.17(d) that could be confused with licensed royalty-bearing products and that do not convey necessary patents rights for the expired 5' Nuclease Process and/or dsDNA-Binding Dye Process
 - No license is conveyed with the purchase of this product under any of US Patents Nos. 5,804,375, 5,994,056, 6,171,785, 5,538,848, 5,723,591, 5,876,930, 6,030,787, and 6,258,569, and corresponding patents outside the United States, or any other patents or patent applications, relating to the 5' Nuclease and dsDNA-Binding Dye Processes. For further information contact Outlicensing at Thermo Fisher Scientific, 5791 Van Allen Way, Carlsbad California 92008; outlicensing@lifetech.com.
- RNA interference
 - RNA interference is licensed under European Patents 1144623, 1214945 and foreign equivalents from Alnylam Pharmaceuticals, Inc., Cambridge, USA and is provided only for use in academic and commercial research (excluding the evaluation or characterization of this product as the potential basis for a siRNA-based drug) and not for any other commercial purposes. Information about licenses for commercial use (including discovery and development of siRNA-based drugs) is available from Alnylam Pharmaceuticals, Inc., 300 Third Street, Cambridge, MA 02142, USA.
- Elitech Dyes and Quenchers
 - Use of this product (including but not limited to Yakima Yellow® and Eclipse™ Dark Quencher) is covered by patents owned by ELITechGroup, Inc., including but not limited to the following US patents and corresponding patent claims outside the US: US 6,699,975, US 6,790,945, US 7,662,942, US 6,653,473, US 6,972,339, 7,541,454, 7,671,218, 7,767,834, 8,163,910, 6,727,356. This product, including both Research Grade Product and IVD Grade Product, is sold under a licensing arrangement with ELITechGroup, Inc. for the purchaser's own internal research. Commercial use of this product requires a separate license from Kaneka Eurogentec S.A. Further information on purchasing licenses may be obtained by contacting licensing@eurogentec.com..
- Custom Molecular Beacon Probes
 - The use of Custom Molecular Beacon Probes employs the following patent rights licensed to Kaneka Eurogentec S.A. from The Public Health Research Institute of the City of New York: (a) the claims of U.S. Patent Application Serial No. 08/439,819; (b) the claims of U.S. Patent Application Serial No. 08/990,176; and (c) all patent application claims and all patent claims, U.S. and foreign, that cover improvements in the foregoing claimed probes, kits and assays and whose subjects are inventions and discoveries that are made by the Public Health Research Institute of New York. The use of Custom Molecular Beacon Probes is limited to the field of use comprising the internal use by an end user of a product solely in assays of the end user (or in applications of the end user's customer, if the end user is performing contract research) in scientific research and development.
- LNA' Oligonucleotides
 - LNA' Oligonucleotides produced under license from Exiqon A/S for research use only. Not for resale or for therapeutic use or use in humans.
- Alexa Fluor®, BODIPY®, Cascade Blue®, Marina Blue®, Oregon Green®, Pacific Blue®, Rhodamine Red™ Texas Red' dyes
 - This product is provided under an intellectual property license from Life Technologies Corporation. The transfer of this product is contingent on the buyer using the purchased product solely in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, 5791 Van Allen Way, Carlsbad, CA 92008 USA or outlicensing@lifetech.com.
- Black Hole Quencher'
 - Black Hole Quencher' and BHQ' are registered trademarks of Biosearch Technologies, Inc., Novato, California, U.S.A. (BTI) Patents are currently pending for the BHQ technology and such BHQ technology is licensed by the manufacturer pursuant to an agreement with BTI and these products are sold exclusively for research and development use only. They may not be used for human or veterinary in vitro or clinical diagnostic purposes and they may not be re-sold, distributed or re-packaged. For information on licensing programs to permit use for human or veterinary in vitro or clinical diagnostic purposes, please contact Biosearch at licensing@biosearchtech.com.
- C-5 Propyne and G-clamp modified Oligonucleotides
 - C-5 Propyne and G-clamp modified oligonucleotides are covered by patents or patents pending owned by Isis Pharmaceuticals, Inc. ("Isis"). Purchase of this product includes a limited license to use this product solely for internal research. This license specifically excludes (and you have no right to use this product for): (a) therapeutic or diagnostic applications (including products or services that incorporate this product), (b) any in vivo toxicity/safety study in support of an investigational new drug application (or foreign counterpart), (c) resale (including sale of any products or services that incorporate this product) or (d) gene functionalization activities (including products or services that

incorporate data derived from gene functionalization activities) if such activities have commercial application, any and all of which require a separate license from Isis. Neither this product nor any product created through its use may be used in human clinical trials. In the event you have a separate agreement with Isis regarding this product which explicitly states that the foregoing is not applicable to you, your use of this product will be governed by the terms of such agreement. In no event does the limited license included with the purchase of this product expand or alter the scope of the license granted pursuant to such agreement.

- PNA
 - All custom PNA oligomers sold by Kaneka Eurogentec S.A. are manufactured by Panagene under an exclusive license and to be used for internal research only. And such custom PNA oligomers are not to be resold unless properly licensed to do so by separate license. All custom PNA oligomers are made and sold pursuant to license under one or more of US Patents Nos. US 5,773,571, 6,133,444, US 6,172,226, US 6,395,474, US 6,414,112, US 6,613,873, US 6,710,163 and US 6,713,602, or corresponding patent claims outside the US.
- Probes
 - Use of this product is covered by one or more of the following US patents and corresponding patent claims outside the US: 5,538,848, 5,723,591, 5,876,930, 6,030,787, 6,258,569 and 5,804,375 (claims 1-12 only). The purchase of this product includes a limited, non-transferable immunity from suit under the foregoing patent claims for using only this amount of product for the purchaser's own internal research. Except under separate license rights available from Applied Biosystems, no right under any other patent claim, or to perform commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, or to sublicense, repackage with other products, or resell in any form, is conveyed expressly, by implication, or by estoppel. This product is for research use only. Diagnostic uses under Roche patents require a separate license from Roche. Further information on purchasing licenses may be obtained from the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.
- Dual-labeled probes or products containing dual-labelled probes that are not Licensed Probes
 - No license is conveyed with the purchase of this product under any of US Patents Nos. 5,804,375, 5,994,056, 6,171,785, 5,538,848, 5,723,591, 5,876,930, 6,030,787, and 6,258,569, and corresponding patents outside the United States, or any other patents or patent applications, relating to the 5' Nuclease and dsDNA-Binding Dye Processes. For further information contact the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.
- HiLyte™ dyes
 - Use of this product is covered by one or more of the following US patents owned by AnaSpec, Inc.: USP 7,465,810, USP 7,754,893, USP 7,820,783, USP 8,258,292, and divisionals, continuations, continuations-in-part, reissues, substitutes, and extensions thereof. The purchase of this product conveys to the buyer a limited, non-exclusive, nontransferable right (without the right to sell, repackage, or further sublicense) under these patent rights to use only the amount of product purchased for the purchaser's own internal research. No other license is granted to the buyer whether expressly, by implication, by estoppel or otherwise. In particular, the purchase of this product does not include nor carry any right or license to use, develop, or otherwise exploit this product commercially, and no rights are conveyed to the buyer to use the product or components of the product for any other purposes, including without limitation, provision of services to a third party, generation of commercial databases, or clinical diagnostics or therapeutics. This product is for research use only. For information on purchasing a license to this product for purposes other than research, contact licensing@eurogentec.com.
- QXL' Quenchers
 - Use of this product is covered by one or more of the following US patents owned by AnaSpec, Inc.: USP 7,910,753, USP 8,093,411, and divisionals, continuations, continuations-in-part, reissues, substitutes, and extensions thereof. The purchase of this product conveys to the buyer a limited, non-exclusive, nontransferable right (without the right to sell, repackage, or further sublicense) under these patent rights to use only the amount of product purchased for the purchaser's own internal research. No other license is granted to the buyer whether expressly, by implication, by estoppel or otherwise. In particular, the purchase of this product does not include nor carry any right or license to use, develop, or otherwise exploit this product commercially, and no rights are conveyed to the buyer to use the product or components of the product for any other purposes, including without limitation, provision of services to a third party, generation of commercial databases, or clinical diagnostics or therapeutics. This product is for research use only. For information on purchasing a license to this product for purposes other than research, contact licensing@eurogentec.com.
- Minor Groove Binder and Eclipse™ Dark Quencher
 - Use of this product is covered by patents owned by ELITechGroup, Inc., including but not limited to the following US patents and corresponding patent claims outside the US: US 6,492,346, US 6,884,584, US 11/784,539, US 7,556,923, US 12/172,999, US-2009-0048427-A1, US 6,084,102, US 6,426,408, US 6,727,356, US 6,790,945, US 6,699,975, US 6,486,308 and US 7,794,945. This product, including both Research Grade Product and IVD Grade Product, is sold under a licensing arrangement with ELITechGroup, Inc. for use by purchaser only in human in vitro diagnostics applications to provide information for individual patient management solely using 5' nuclease (Taqman-MGB) technology. Commercial use of this product requires a separate license from Kaneka Eurogentec S.A. Further information on purchasing licenses may be obtained by contacting licensing@eurogentec.com.
- MOE
 - The use of this product is covered by patents or patents pending owned by Ionis Pharmaceuticals, Inc. ("Ionis"). Purchase of this product includes a limited, non-transferable right to use the purchased amount of the product for the purchaser's own internal research and does not include any license to use this product for other purpose like (for instance and not exhaustive): (a) therapeutic or diagnostic applications, (b) any in vivo toxicity/safety study in support of an investigational new drug application (or foreign counterpart), (c) resale (including sale of any products or services that incorporate this product or oligonucleotides) or (d) gene functionalization activities (including products or services that incorporate data derived from gene functionalization activities). Any and all of such application of the oligonucleotides may require a separate license from Ionis.

Alexa Fluor® is a registered trademark of Molecular Probes, Inc. (a wholly owned subsidiary of Invitrogen Corp.)
 BHQ-0®, BHQ-1®, BHQ-2®, BHQ-3® are registered trademarks of Biosearch Technologies, Inc.
 Black Hole Quencher® is a registered trademark of Biosearch Technologies, Inc.
 CAL Fluor™ is trademark of Biosearch Technologies, Inc.
 BODIPY® is a registered trademark of Molecular Probes, Inc. (a wholly owned subsidiary of Invitrogen Corp.)
 Cascade Blue® is a registered trademark of Molecular Probes, Inc. (a wholly owned subsidiary of Invitrogen Corp.)
 Cy® is a registered trademark of GE Healthcare Bio-Sciences Corp.
 Cy®5 Direct is a registered trademark of GE Healthcare Bio-Sciences Corp.
 Dragonfly Orange™ is a trademark of Kaneka Eurogentec S.A.
 Eclipse® is a registered trademark of Epoch Biosciences, Inc.
 Excel™ is a trademark of Microsoft Corp.
 EGT™ is a trademark of Kaneka Eurogentec S.A.
 Eurogentec® is a registered trademark of Kaneka Eurogentec S.A.
 HiLyte™ is a trademark of Anaspec, Inc.
 LNA® is a registered trademark of Exiqon A/S.
 Marina Blue® is a registered trademark of Molecular Probes, Inc. (a wholly owned subsidiary of Invitrogen Corp.)
 MGB® is a registered trademark of Epoch Biosciences, Inc.
 MGB Eclipse® is a registered trademark of Epoch Biosciences, Inc.
 Oli&GO™ is a trademark of Kaneka Eurogentec S.A.
 Oregon Green® is a registered trademark of Molecular Probes, Inc. (a wholly owned subsidiary of Invitrogen Corp.)
 Pacific Blue™ is a trademark of Molecular Probes, Inc. (a wholly owned subsidiary of Invitrogen Corp.)
 Plexor™ is a trademark of Promega Corporation
 Quasar™ is trademark of Biosearch Technologies, Inc.
 QXL® is a registered trademark of Anaspec, Inc.
 Rhodamine Red™ and Rhodamine Green™ are trademarks of Molecular Probes, Inc. (a wholly owned subsidiary of Invitrogen Corp.)
 ROX™ is a trademark of Applera Corp.
 RP-Cartridge-Gold™ is a trademark of Kaneka Eurogentec S.A.
 Taqman® is a registered trademark of Applied Biosystems, Inc.
 Texas Red® is a registered trademark of Molecular Probes, Inc. (a wholly owned subsidiary of Invitrogen Corp.)
 Track™ is a trademark of Kaneka Eurogentec S.A.
 Unique™ is a trademark of Kaneka Eurogentec S.A.
 UltraPureGold™ is a trademark of Kaneka Eurogentec S.A.
 Vic® is a registered trademark of Applera Corp.
 Yakima Yellow® is a registered trademark of Epoch Biosciences, Inc.

www.eurogentec.com

Oligos

Contact:

For more information,
please contact your Sales
representative or contact us at:
oligo@eurogentec.com

Kaneka Eurogentec S.A.

rue Bois Saint-Jean, 5
4102 Seraing
Belgium

Toll free number: **00 800 666 00 123**

Eurogentec 2018

Follow us on



Life Science Research



Oligos




PCR/qPCR



Genes/DNA



Antibodies



Peptides



Proteins