

# pipetman®

L



## VARIABLE VOLUME MECHANICAL PIPETTES



**PIPETMAN® L multichannel models now cover a volume range of 0.5 µL to 1200 µL.**

### I Ultimate Performance

Built for traceability and comfort, PIPETMAN L is designed to increase daily productivity, giving you the ability to pipette longer without fatigue.

### I Secure Pipetting

Reduce pipetting errors and avoid accidental volume changes during pipetting cycles with a volume-locking mechanism for improved accuracy and precision.

### I Optimized Comfort

Lightweight with light pipetting forces. Right-handed or left-handed? No problem. Adjust the location of the tip ejector to the most comfortable position for you.





## Classic metal or plastic tip ejector

Choose between a stainless steel tip ejector for more robustness or a plastic tip ejector to avoid corrosion.



## Choose your model

Choose your channel number and volume range. Each ensures consistent volume levels, allowing you to aspirate and dispense your samples with confidence.

**NEW**

### PIPETMAN L Single Channel — Variable Volume Models

Model	PIPETMAN® DIAMOND Tips		Part Number	Volume (μL)
P2L	D10 DL10	DF10ST DFL10ST	FA10001P	0.2 – 2
			FA10001M	
P10L*	D10 DL10	DF10ST DFL10T	FA10002P	1 – 10
			FA10002M	
P20L	D200	DF30ST	FA10003P	2 – 20
			FA10003M	
P100L	D200	DF100ST	FA10004P	10 – 100
			FA10004M	
P200L	D200 D300	DF200ST DF300ST	FA10005P	20 – 200
			FA10005M	
P1000L	D1000 D1200	DF1000ST DF1200ST	FA10006P	100 – 1000
			FA10006M	
P5000L	D5000		FA10007**	500 – 5000
P10mLL	D10mL		FA10008**	1000 – 10000

### PIPETMAN L Multichannel

<b>P8x10L</b>	D10 DL10	DF10ST DFL10ST	<b>FA10013</b>	0.5 – 10
<b>P12x10L</b>			<b>FA10014</b>	
<b>P8x20L</b>	DL10 D200	DFL10ST DF30ST	<b>FA10009</b>	2 – 20
<b>P12x20L</b>			<b>FA10010</b>	
<b>P8x200L</b>	D200 D300	DF200ST DF300ST	<b>FA10011</b>	20 – 200
<b>P12x200L</b>			<b>FA10012</b>	
<b>P8x300L</b>	D200 D300	DF200ST DF300ST	<b>FA10015</b>	20 – 300
<b>P12x300L</b>			<b>FA10016</b>	
<b>P8x1200L</b>	D1200	DF1200ST	<b>FA10039</b>	100 – 1200
<b>P12x1200L</b>			<b>FA10040</b>	

\*P10L model can be used up to 0.5 μL

For a pipette with a plastic tip ejector, the part number ends with the letter P.

For a pipette with a stainless steel tip ejector, the part number ends with the letter M.

PIPETMAN L specifications are downloadable from [www.gilson.com](http://www.gilson.com)

PIPETMAN L single channel models are available in kits.

For even more reproducibility, PIPETMAN L is also available in a fixed version.

[gilson.com/contactus](http://gilson.com/contactus)

**ISO 9001**  
CERTIFIED

LT800565/B | ©2020 Gilson S.A.S.  
All rights reserved. All trademarks  
are the property of Gilson, Inc.  
and its subsidiaries.

**GILSON®**